

Imaging·Illumination·Spectroscopy·Microscopy

Components and Tools 2011 Catalog



TET-EPI
ILLUMINATOR

WORLDWIDE SHIPMENT · TIMELY TECHNICAL SUPPORT



光技術をサポートする

株式会社オプトサイエンス

<http://www.optoscience.com>

東京本社 〒160-0014 東京都新宿区内藤町1番地 内藤町ビルディング
TEL:03(3356)1064 FAX:03(3356)3466 E-mail:info@optoscience.com
大阪支店 〒532-0011 大阪市淀川区西中島7-7-2 新大阪ビル西館
TEL:06(6305)2064 FAX:06(6305)1030 E-mail:osk@optoscience.com
名古屋営業所 〒450-0002 名古屋市中村区名駅2-37-21 東海ソフトビル
TEL:052(569)6064 FAX:052(569)8064 E-mail:ngo@optoscience.com

ORDER ONLINE

www.mightex.com
www.mightexsystems.com

 **Mightex**
Simply Brighter

ORDERING INFORMATION

Five Easy Ways to Order -

1. **Website** — order directly online at www.mightex.com or www.mightexsystems.com.



Place your order online, or call
(USA) +1-925-218-1885
Worldwide Shipment



Order online!
Ship today!
1st-class support!

2. **E-mail** — sales@mightex.com or sales@mightexsystems.com.
3. **Telephone** — (USA) 1-925-218-1885 or (Canada) 1-416-840-4991
4. **Fax** — 1-416-840-6541.
5. **Mail** — Mightex Systems
2343 Brimley Road
P.O. Box 868
Toronto, Ontario M1S 3L6
CANADA

Payment Terms -

1. **US and Canadian customers** — eligible of NET-30 terms, upon credit approval;
2. **International customers** — prepayment (e.g. credit card, wire transfer or money order) is generally required;
3. **Credit cards and other payment methods** — VISA, MasterCard, Discover, American Express, PayPal and eCheck.

Shipping -

All shipments are made from California, unless otherwise stated. A street address is **required**. If you have a post office box **only**, we will ship via U.S. Postal Service.

Shipping cost estimate? — Please refer to **Quotation** below.

Quotation -

You may obtain an official quotation (including shipping cost) directly on www.mightex.com or www.mightexsystems.com. Simply add products to shopping cart, and then click "Get Quote". You would need to provide your exact shipping address, and our eCommerce system would be able to calculate the shipping cost.

Prices -

All product prices are posted directly on www.mightex.com or www.mightexsystems.com, and they are subject to change without notice. We also reserve the right to correct misprints, and will notify you before charging you more than the price listed on Mightex's website.

Specifications -

Listed specifications are accurate as of the publication date, and they may change without prior notice due to product improvements and design changes.

Product Returns -

If you need to return a product, please contact us at sales@mightex.com for an authorization number (i.e. RMA#) to aid in tracking returned goods. For more details, please refer to Mightex's **Cancellation and Return Policy** below.

Cancellation and Return Policy -

Please refer to: <http://www.mightexsystems.com/support/cancellationsandreturnpolicy.html>.

Warranty Terms and Conditions -

Please refer to: <http://www.mightexsystems.com/support/warrantytermsandconditions.html>.

■ Industrial Cameras

CMOS Camera.....	06
CCD Cameras	17
Line Scan Cameras.....	25

■ LED Controllers

Manual Controlled	30
Computer Controlled	32
Manual & Computer (“Dual Mode”) Controlled	37

■ Light Sources

Fiber-Coupled LED’s	40
Light-Guide-Coupled LED’s.....	41
LED Collimator Sources	42
Low-Cost LED Spot Lights	43
Precision LED Spot Lights	45
Spectrum Synthesizing Sources (“Cubic-S”)	46

■ Microscope Tools

Microscopy Light Sources	48
Low-Cost LED Based Epi-Fluorescence Microscopy Solution	49
Fluorescence Illuminators for stereo-Microscopes	50
Filter Reader Spectrophotometer	51

■ Spectrometers

High-Resolution High-Stability CCD Spectrometers.....	52
---	----




Table of Content

USB2.0 Area Cameras

Unbuffered USB2.0 CMOS Cameras (M-Series)	6
Buffered USB2.0 CMOS Cameras (B-Series)	8
Buffered USB2.0 CCD Cameras (C-Series)	17
Windowless USB2.0 CMOS Cameras (W-Series).....	23

USB2.0 Line Cameras

3648 Pixel 16-bit CCD Line Camera with External Trigger	25
High-Speed 2048-Pixel 12-bit CCD Line Camera with External Trigger	27

Unbuffered USB2.0 Monochrome 1.3MP CMOS Cameras (M-Series)

Mightex USB 2.0 cameras are designed for a wide variety of applications (such as industrial inspections, digital microscopy and medical imaging), which require good quality cameras that are easy to use and cost-effective. Monochrome cameras often exhibit 20% higher spatial resolution than their color counterparts because no pixel interpolation is necessary. Without a Bayer color filter on the sensor, monochrome cameras are also more sensitive than color sensors especially in near IR and UV regions. Frame rate can be as high as 22 fps at full resolution and up to 600 fps using ROI mode. In addition, a user-friendly GUI based application software and an SDK are provided for custom software development. The cameras have 4-pin GPIOs and optional built-in LED drivers. A DirectShow driver and a TWAIN driver are also provided to easily link the cameras with users' applications.



PERFORMANCE SPECIFICATIONS

Parameters	MCN-B013-U	GLN-B013-U	MCE-B013-U	MLE-B013-U	Unit
Board-level/enclosed	Board-level	Board-level	Enclosed	Enclosed	
Number of GPIOs	4	4	4	No	
Built-in LED Drivers	No	Yes	No	Yes	
Resolution	1280x1024 monochrome				
CMOS Chip	½" (5:4) Aptina MT9M001, Rolling Shutter				
Pixel Size	5.2 x 5.2				µm
Active Imager Size	6.66 x 5.32				mm
Scanning System	Progressive				
Dynamic Range	68				dB
Sensor SNR	45				dB
Gray Level	8				bit
Responsivity	2.1				V/lux -sec
Frame Rates (@48MHz Clock)	22 @1280x1024, 32 @1024x768, 50 @800x600, 70 @640x480 180 @320x240, 300 @160x120, 450 @64x64, 600 @32x32				fps
Sub Resolutions	1024x768, 800x600, 640x480, 320x240, 160x120, 64x64, 32x32				
Shutter Speed (Exposure time)	0.04~750				ms
Hardware Gains	12				dB
Trigger Mode	With external trigger				
Lens mount	C- mount or CS-mount (M12.5-mount or custom-defined lens mount supported)				
Built-in Filters	IR-cut (factory standard), No filter, or IR-pass				
Power consumption	< 1.0 (excluding LED drivers, if applicable)				W
Number of LED Driver Channels	N.A.	4	N.A.	4	
LED Driver Max. Output Voltage	N.A.	5	N.A.	5	V
LED Driver Max. Output Current (total)	N.A.	250	N.A.	250	mA
Dimension	51 x 51 x 29 (CS-mount) 51 x 51 x 34 (C-mount)		58 x 58 x 34 (CS-mount) 58 x 58 x 39 (C-mount)		mm
Weight (excluding lens)	48		150		g

Unbuffered USB2.0 Color 3MP CMOS Cameras (M-Series)

Mightex USB 2.0 cameras are designed for a wide variety of applications (such as industrial inspections, digital microscopy and medical imaging), which require good quality cameras that are easy to use and cost-effective. With a USB2.0 interface and powerful PC software, the camera delivers excellent quality images, and the frame rate can be as high as 7 fps at full resolution and up to 600 fps using ROI mode. In addition, a user-friendly GUI based application software and a SDK are provided for custom software development. The cameras have 4-pin GPIOs and optional built-in LED drivers. A DirectShow driver and a TWAIN driver are also provided to easily link the cameras with users' applications.



PERFORMANCE SPECIFICATIONS

Parameters	MCN-C030-U	GLN-C030-U	MCE-C030-U	MLE-C030-U	Unit
Board-level/enclosed	Board-level	Board-level	Enclosed	Enclosed	
Number of GPIOs	4	4	4	No	
Built-in LED Drivers	No	Yes	No	Yes	
Resolution	2048x1536 color				
CMOS Chip	1/2" Aptina MT9T001, rolling shutter				
Pixel Size	3.2 x 3.2				μm
Scanning System	Progressive				
Dynamic Range	61				dB
Sensor SNR	43				dB
Color (RGB)	8				bit
Frame Rates (@48MHz Clock)	7 @2048x1536, 11 @1600x1200, 16 @1280x1024, 26 @1024x768, 35 @800x600, 50 @640x480, 80 @320x240, 130 @160x120, 175 @64x64				fps
Sub Resolutions	1600x1200, 1280x1024, 1024x768, 800x600, 640x480, 320x240, 160x120, 64x64				
Shutter Speed (Exposure time)	0.04~750				ms
Hardware Gains	12				dB
Trigger Mode	With external trigger				
Lens mount	C- mount or CS-mount (M12.5-mount or custom-defined lens mount supported)				
Built-in Filters	IR-cut (factory standard), IR-pass or no filter				
Power consumption	< 1.0 (excluding LED drivers, if applicable)				W
Number of LED Driver Channels	N.A.	4	N.A.	4	
LED Driver Maximum Output Voltage	N.A.	5	N.A.	5	V
LED Driver Maximum Output Current (total)	N.A.	250	N.A.	250	mA
Dimension	51 x 51 x 29 (CS-mount) 51 x 51 x 34 (C-mount)		58 x 58 x 34 (CS-mount) 58 x 58 x 39 (C-mount)		mm
Weight (excluding lens)	48		150		g

Buffered USB2.0 Mono 752x480 CMOS Cameras 8-bit (B-Series)

Mightex USB2.0 cameras with frame buffers are optimized for machine-vision applications, and they can be also used for a wide variety of applications (such as industrial inspections, digital microscopy and medical imaging), which require good quality cameras that are easy to use and cost-effective. These cameras have built-in frame buffers, external trigger-in, strobe-out, and a powerful camera engine that supports multiple cameras. Monochrome cameras often exhibit 20% higher spatial resolution than their color counterparts because no pixel interpolation is necessary. Since there is no *Bayer color filter* on the sensor, monochrome cameras are more sensitive than color sensors, especially in near IR and UV regions. Frame rate can be as high as 60 fps in full resolution and up to 600 fps using ROI mode. A user-friendly GUI based application software and an SDK are provided for custom software development. In addition, a DirectShow driver and a TWAIN driver are available to easily link the cameras with users' applications. A USB command set protocol and a Linux driver are also provided for non-Windows based applications.



PERFORMANCE SPECIFICATIONS

Parameters	BCN-BG04-U	BCE-BG04-U	Unit
Board-level/enclosed	Board-level	Enclosed	
Number of GPIOs	4	4	
Resolution	752x480 Monochrome		
CMOS Chip	1/3" Aptina MT9V032, global shutter		
Pixel Size	6.0x6.0		μm
Scanning System	Progressive		
Dynamic Range	>55		dB
Gray Level	8		bit
Responsivity	4.8		V/lux-sec
On-board Memory ("Frame Buffers")	32		MB
Frame Rates(@26MHz Clock)	60 @752x480, 65 @640x480, 130 @320x240 220 @160x120, 310 @64 x 64		fps
Sub Resolutions	Support arbitrary ROI (Nx, Ny), with Nx and Ny multiples of 4		
Shutter Speed (Exposure time)	0.05 ~ 750		ms
Hardware Gains	1x ~ 4x		
Trigger Mode	With external trigger		
Strobe Out	Yes		
Lens mount	C- mount or CS-mount (M12.5-mount or custom-defined lens mount supported)		
Built-in Filters	IR-cut (factory standard), or IR-pass, or no filter		
Power consumption	< 1.8		W
Dimension	51 x 51 x 29 (CS-mount) 51 x 51 x 34 (C-mount)	58 x 58 x 34 (CS-mount) 58 x 58 x 39 (C-mount)	mm
Weight (excluding lens)	80	150	g

Buffered USB2.0 Color 752x480 CMOS Cameras (B-Series)

Mightex USB 2.0 cameras are designed for a wide variety of applications (such as industrial inspections, digital microscopy and medical imaging) which require good quality cameras that are easy to use and cost-effective. With a USB2.0 interface and powerful PC software, the camera delivers excellent quality images, and the frame rate can be as high as 60 fps in full resolution and up to 600 fps using ROI mode. In addition, a user-friendly GUI based application software and a SDK are provided for custom software development. A DirectShow driver and a TWAIN driver are available to easily link the cameras with users' applications. For non-Windows based applications, a USB command set protocol and a Linux driver are also provided.



PERFORMANCE SPECIFICATIONS

Parameters	BCN-CG04-U	BCE-CG04-U	Unit
Board-level/enclosed	Board-level	Enclosed	
Number of GPIOs	4	4	
Resolution	752x480 color		
CMOS Chip	1/3" Aptina MT9V032, global shutter		
Pixel Size	6.0x6.0		μm
Scanning System	Progressive		
Dynamic Range	>55		dB
Gray Level	8		bit
Responsivity	4.8		V/lux-sec
On-board Memory ("Frame Buffers")	32ddccvd		MB
Frame Rates (@26MHz Clock)	60 @752x480, 65 @640x480 130 @320x240, 220 @160x120, 310 @64x64		fps
Sub Resolutions	Support arbitrary ROI (Nx, Ny), with Nx and Ny multiples of 4		
Shutter Speed (Exposure time)	0.05 ~ 750		ms
Hardware Gains	1x ~ 4x		
Trigger Mode	With external trigger		
Strobe Out	Yes		
Lens mount	C- mount or CS-mount (M12.5-mount or custom-defined lens mount supported)		
Built-in Filters	IR-cut (factory standard), or IR-pass, or no filter		
Power consumption	< 1.8		W
Dimension	51 x 51 x 29 (CS-mount) 51 x 51 x 34 (C-mount)	58 x 58 x 34 (CS-mount) 58 x 58 x 39 (C-mount)	mm
Weight (excluding lens)	80	150	g

Buffered USB2.0 Monochrome 752x480 CMOS Cameras 10-bit (B-Series)

Mightex USB2.0 cameras with frame buffers are optimized for machine-vision applications, and they can be also used for a wide variety of applications (such as industrial inspections, digital microscopy and medical imaging), which require good quality cameras that are easy to use and cost-effective. These cameras have built-in frame buffers, external trigger-in, strobe-out, and a powerful camera engine that supports multiple cameras. Monochrome cameras often exhibit 20% higher spatial resolution than their color counterparts because no pixel interpolation is necessary. Since there is no *Bayer color filter* on the sensor, monochrome cameras are more sensitive than color sensors, especially in near IR and UV regions. Frame rate can be as high as 60 fps in full resolution and up to 600 fps using ROI mode. A user-friendly GUI based application software and an SDK are provided for custom software development. In addition, a Direct-Show driver and a TWAIN driver are available to easily link the cameras with users' applications. A USB command set protocol and a Linux driver are also provided for non-Windows based applications.



PERFORMANCE SPECIFICATIONS

Parameters	BTN-BG04-U		BTE-BG04-U	Unit
Board-level/enclosed	Board-level		Enclosed	
Number of GPIOs	4		4	
Resolution	752x480 Monochrome			
CMOS Chip	1/3" Aptina MT9V032, global shutter			
Pixel Size	6.0 x 6.0			μm
Scanning System	Progressive			
Dynamic Range	>55			dB
Gray Level	10			bit
Responsivity	4.8			V/lux-sec
On-board Memory ("Frame Buffers")	32			MB
Frame Rates (@26MHz Clock)	8 bit operation mode	60 @752x480, 65 @640x480, 130 @320x240 220 @160x120, 310 @64 x 64		fps
	10 bit operation mode	36 @752x480, 44 @640x480, 130 @320x240 220 @160x120, 310 @64 x 64		
Sub Resolutions	Support arbitrary ROI (Nx, Ny), with Nx and Ny multiples of 4			
Shutter Speed (Exposure time)	0.05 ~ 750			ms
Hardware Gains	1x ~ 4x			
Trigger Mode	With external trigger			
Strobe Out	Yes			
Lens mount	C– mount or CS-mount (M12.5-mount or custom-defined lens mount supported)			
Built-in Filters	IR-cut (factory standard), or IR-pass, or no filter			
Power consumption	< 1.8			W
Dimension	51 x 51 x 29 (CS-mount) 51 x 51 x 34 (C-mount)	58 x 58 x 34 (CS-mount) 58 x 58 x 39 (C-mount)		mm
Weight (excluding lens)	80		150	g

Buffered USB2.0 Monochrome 1.3MP CMOS Cameras 8-bit (B-Series)

Mightex USB2.0 cameras with frame buffers are optimized for machine-vision applications, and they can be also used for a wide variety of applications (such as industrial inspections, digital microscopy and medical imaging), which require good quality cameras that are easy to use and cost-effective. These cameras have built-in frame buffers, external trigger-in, strobe-out, and a powerful camera engine that supports multiple cameras. Monochrome cameras often exhibit 20% higher spatial resolution than their color counterparts because no pixel interpolation is necessary. Since there is no *Bayer color filter* on the sensor, monochrome cameras are more sensitive than color sensors, especially in near IR and UV regions.

Frame rate can be as high as 25 fps in full resolution and up to 300 fps using ROI mode. A user-friendly GUI based application software and an SDK are provided for custom software development. In addition, a DirectShow driver and a TWAIN driver are available to easily link the cameras with users' applications. A USB command set protocol and a Linux driver are also provided for non-Windows based applications.



PERFORMANCE SPECIFICATIONS

Parameters	BCN-B013-U	BCE-B013-U	Unit
Board-level/enclosed	Board-level	Enclosed	
Number of GPIOs	4	4	
Resolution	1280x1024 monochrome		
CMOS Chip	½" (5:4) Aptina MT9M001, rolling shutter		
Pixel Size	5.2 x 5.2		μm
Active Imager Size	6.66 x 5.32		mm
Scanning System	Progressive		
Dynamic Range	68		dB
Sensor SNR	45		dB
Gray Level	8		bit
Responsivity	2.1		V/lux-sec
On-board Memory ("Frame Buffers")	32		MB
Frame Rates (@48MHz Clock)	25 @1280x1024, 35 @1024x768, 50 @800x600, 55 @752x480, 65 @640x480, 140 @320x240		fps
Sub Resolutions	Support arbitrary ROI (Nx, Ny), with Nx and Ny multiples of 4		
Shutter Speed (Exposure time)	0.05 ~ 750		ms
Hardware Gains	0.125x ~ 8x		
Trigger Mode	With external trigger		
Strobe Out	Yes		
Lens mount	C- mount or CS-mount (M12.5-mount or custom-defined lens mount supported)		
Built-in Filters	IR-cut (factory standard), IR-pass, or no filter		
Power consumption	< 1.8		W
Dimension	51 x 51 x 29 (CS-mount) 51 x 51 x 34 (C-mount)	58 x 58 x 34 (CS-mount) 58 x 58 x 39 (C-mount)	mm
Weight (excluding lens)	80	150	g

Buffered USB2.0 Monochrome 1.3MP CMOS Cameras 10-bit (B-Series)

Mightex USB2.0 cameras with frame buffers are optimized for machine-vision applications, and they can be also used for a wide variety of applications (such as industrial inspections, digital microscopy and medical imaging), which require good quality cameras that are easy to use and cost-effective. These cameras have built-in frame buffers, external trigger-in, strobe-out, and a powerful camera engine that supports multiple cameras. Monochrome cameras often exhibit 20% higher spatial resolution than their color counterparts because no pixel interpolation is necessary. Since there is no *Bayer color filter* on the sensor, monochrome cameras are more sensitive than color sensors, especially in near IR and UV regions. Frame rate can be as high as 12 fps in full resolution and up to 150 fps using ROI mode. A user-friendly GUI based application software and an SDK are provided for custom software development. In addition, a DirectShow driver and a TWAIN driver are available to easily link the cameras with users' applications. A USB command set protocol and a Linux driver are also provided for non-Windows based applications.



PERFORMANCE SPECIFICATIONS

Parameters	BTN-B013-U	BTE-B013-U	Unit
Board-level/enclosed	Board-level	Enclosed	
Number of GPIOs	4	4	
Resolution	1280x1024 Monochrome		
CMOS Chip	½" (5:4) Aptina MT9M001, Rolling Shutter		
Pixel Size	5.2 x 5.2		µm
Active Imager Size	6.66 x 5.32		mm
Scanning System	Progressive		
Dynamic Range	68		dB
Sensor SNR	45		dB
Gray Level	10		bit
Responsivity	2.1		V/lux-sec
On-board Memory ("Frame Buffers")	32		MB
Frame Rates (@48MHz Clock)	12 @1280x1024, 18 @1024x768, 25 @800x600 28 @752x480, 33 @640x480, 70 @320x240		fps
Sub Resolutions	Support arbitrary ROI (Nx, Ny), with Nx and Ny multiples of 4		
Shutter Speed (Exposure time)	0.05 ~ 750		ms
Hardware Gains	0.125x ~ 8x		
Trigger Mode	With external trigger		
Strobe Out	Yes		
Lens mount	C- mount or CS-mount (M12.5-mount or custom-defined lens mount supported)		
Built-in Filters	IR-cut (factory standard), IR-pass, or no filter		
Power consumption	< 1.8		W
Dimension	51 x 51 x 29 (CS-mount) 51 x 51 x 34 (C-mount)	58 x 58 x 34 (CS-mount) 58 x 58 x 39 (C-mount)	mm
Weight (excluding lens)	80	150	g

Buffered USB2.0 Color 3MP CMOS Cameras (B-Series)

Mightex USB 2.0 cameras are designed for a wide variety of applications (such as industrial inspections, digital microscopy and medical imaging), which require good quality cameras that are easy to use and cost-effective. With a USB2.0 interface and powerful PC software, the camera delivers excellent quality images, and the frame rate can be as high as 8 fps in full resolution and up to 300 fps using ROI mode. In addition, a user-friendly GUI based application software and a SDK are provided for custom software development. A DirectShow driver and a TWAIN driver are available to easily link the cameras with users' applications. For non-Windows based applications, a USB command set protocol and a Linux driver are also provided.



PERFORMANCE SPECIFICATIONS

Parameters	BCN-C030-U	BCE-C030-U	Unit
Board-level/enclosed	Board-level	Enclosed	
Number of GPIOs	4	4	
Resolution	2048 x 1536 color		
CMOS Chip	1/2" (5:4) Aptina MT9T001, rolling shutter		
Pixel Size	3.2 x 3.2		μm
Scanning System	Progressive		
Dynamic Range	61		dB
Sensor SNR	43		dB
Gray Level	8		bit
Responsivity	2.1		V/lux-sec
On-board Memory ("Frame Buffers")	32		MB
Frame Rates (@48MHz Clock)	8 @2048x1536, 14 @1600x1200, 25 @1280x1024 35 @1024x768, 50 @800x600, 55 @752x480 65 @640x480, 140 @320x240		fps
Sub Resolutions	Support arbitrary ROI (Nx, Ny), with Nx and Ny multiples of 4		
Shutter Speed (Exposure time)	0.05 ~ 750		ms
Hardware Gains	0.125x ~ 8x		
Trigger Mode	With external trigger		
Strobe Out	Yes		
Lens mount	C- mount or CS-mount (M12.5-mount or custom-defined lens mount supported)		
Built-in Filters	IR-cut (factory standard), or IR-pass, or no filter		
Power consumption	< 1.8		W
Dimension	51 x 51 x 29 (CS-mount) 51 x 51 x 34 (C-mount)	58 x 58 x 34 (CS-mount) 58 x 58 x 39 (C-mount)	mm
Weight (excluding lens)	80	150	g

Buffered USB2.0 Monochrome 5MP CMOS Cameras 8-bit (B-Series)

Mightex USB2.0 cameras with frame buffers are optimized for machine-vision applications, and they can be also used for a wide variety of applications (such as industrial inspections, digital microscopy and medical imaging) which require good quality cameras that are easy to use and cost-effective. These cameras have built-in frame buffers, external trigger-in, strobe-out, and a powerful camera engine that supports multiple cameras. Monochrome cameras often exhibit 20% higher spatial resolution than their color counterparts because no pixel interpolation is necessary. Since there is no *Bayer color filter* on the sensor, monochrome cameras are more sensitive than color sensors, especially in near IR and UV regions. Frame rate can be as high as 6 fps in full resolution and up to 150 fps using ROI mode. A user-friendly GUI based application software and an SDK are provided for custom software development. In addition, a DirectShow driver and a TWAIN driver are available to easily link the cameras with users' applications. A USB command set protocol and a Linux driver are also provided for non-Windows based applications.



PERFORMANCE SPECIFICATIONS

Parameters	BCN-B050-U	BCE-B050-U	Unit
Board-level/enclosed	Board-level	Enclosed	
Number of GPIOs	4	4	
Resolution	2592 x 1944 Monochrome		
CMOS Chip	Aptina MT9P031, rolling shutter		
Pixel Size	2.2 x 2.2		μm
Scanning System	Progressive		
Dynamic Range	70		dB
Sensor SNR	38		dB
Gray Level	8		bit
Responsivity	1.4		V/lux-sec
Frame Rates (@48MHz Clock)	6 @2592x1944, 9 @2048x1536, 13 @1600x1200, 18 @1280x1024, 28 @1024x768, 40 @800x600, 50 @752x480, 56 @640x480, and 100 @320x240		fps
Sub Resolutions	Support arbitrary ROI (Nx, Ny), with Nx and Ny multiples of 4		
Shutter Speed (Exposure time)	0.05 ~ 750		ms
Hardware Gains	1x ~ 16x		
Trigger Mode	With external trigger		
Strobe Out	Yes		
Lens mount	C– mount and CS-mount (M12.5-mount or custom-defined lens mount supported)		
Built-in Filters	IR-cut (factory standard), or IR-pass, or no filter		
Power consumption	< 1.8		W
Dimension	51 x 51 x 29 (CS-mount) 51 x 51 x 34 (C-mount)	58 x 58 x 34 (CS-mount) 58 x 58 x 39 (C-mount)	mm
Weight (excluding lens)	80	150	g

Buffered USB2.0 Color 5MP CMOS Cameras (B-Series)

Mightex USB 2.0 cameras are designed for a wide variety of applications (such as industrial inspections, digital microscopy and medical imaging), which require good quality cameras that are easy to use and cost-effective. With a USB2.0 interface and powerful PC software, the camera delivers excellent quality images, and the frame rate can be as high as 6 fps in full resolution and up to 600 fps using ROI mode. In addition, a user-friendly GUI based application software and a SDK are provided for custom software development. A DirectShow driver and a TWAIN driver are available to easily link the cameras with users' applications. For non-Windows based applications, a USB command set protocol and a Linux driver are also provided.



PERFORMANCE SPECIFICATIONS

Parameters	BCN-C050-U	BCE-C050-U	Unit
Board-level/enclosed	Board-level	Enclosed	
Number of GPIOs	4	4	
Resolution	2592 x 1944 color		
CMOS Chip	Aptina MT9P031, rolling shutter		
Pixel Size	2.2 x 2.2		μm
Scanning System	Progressive		
Dynamic Range	70		dB
Sensor SNR	38		dB
Gray Level	8		bit
Responsivity	1.4		V/lux-sec
Frame Rates (@48MHz Clock)	6 @2592x1944, 9 @2048x1536, 13 @1600x1200, 18 @1280x1024, 28 @1024x768, 40 @800x600, 50 @752x480, 56 @640x480, 100 @320x240		fps
Sub Resolutions	Support arbitrary ROI (Nx, Ny), with Nx and Ny multiples of 4		
Shutter Speed (Exposure time)	0.05 ~ 750		ms
Hardware Gains	1x ~ 16x		
Trigger Mode	With external trigger		
Strobe Out	Yes		
Lens mount	C- mount or CS-mount (M12.5-mount or custom-defined lens mount supported)		
Built-in Filters	IR-cut (factory standard), or IR-pass, or no filter		
Power consumption	< 1.8		W
Dimension	51 x 51 x 29 (CS-mount) 51 x 51 x 34 (C-mount)	58 x 58 x 34 (CS-mount) 58 x 58 x 39 (C-mount)	mm
Weight (excluding lens)	80	150	g

Buffered USB2.0 Monochrome 5MP CMOS Cameras 12-bit (B-Series)

Mightex USB2.0 cameras with frame buffers are optimized for machine-vision applications, and they can be also used for a wide variety of applications (such as industrial inspections, digital microscopy and medical imaging) which require good quality cameras that are easy to use and cost-effective. These cameras have built-in frame buffers, external trigger-in, strobe-out, and a powerful camera engine that supports multiple cameras. Monochrome cameras often exhibit 20% higher spatial resolution than their color counterparts because no pixel interpolation is necessary. Since there is no *Bayer color filter* on the sensor, monochrome cameras are more sensitive than color sensors, especially in near IR and UV regions. Frame rate can be as high as 6 fps in full resolution and up to 150 fps using ROI mode. A user-friendly GUI based application software and an SDK are provided for custom software development. In addition, a DirectShow driver and a TWAIN driver are available to easily link the cameras with users' applications. A USB command set protocol and a Linux driver are also provided for non-Windows based applications.




PERFORMANCE SPECIFICATIONS

Parameters	BTN-B050-U	BTE-B050-U	Unit
Board-level/enclosed	Board-level	Enclosed	
Number of GPIOs	4	4	
Resolution	2592 x 1944 Monochrome		
Bit	8/12 (switchable)		bits
CMOS Chip	Aptina MT9P031, rolling shutter		
Pixel Size	2.2 x 2.2		μm
Scanning System	Progressive		
Dynamic Range	70		dB
Sensor SNR	38		dB
Gray Level	8		bit
Responsivity	1.4		V/lux-sec
Frame Rates (@48MHz Clock)	3 @2592x1944, 5 @2048x1536, 7 @1600x1200, 12 @1280x1024, 21 @1024x768, 33 @800x600, 43 @752x480, 51 @640x480, and 100 @320x240		fps
Sub Resolutions	Support arbitrary ROI (Nx, Ny), with Nx and Ny multiples of 4		
Shutter Speed (Exposure time)	0.05 ~ 750		ms
Hardware Gains	1x ~ 16x		
Trigger Mode	With external trigger		
Strobe Out	Yes		
Lens mount	C- mount and CS-mount (M12.5-mount or custom-defined lens mount supported)		
Built-in Filters	IR-cut (factory standard), or IR-pass, or no filter		
Power consumption	< 1.8		W
Dimension	51 x 51 x 29 (CS-mount) 51 x 51 x 34 (C-mount)	58 x 58 x 34 (CS-mount) 58 x 58 x 39 (C-mount)	mm
Weight (excluding lens)	80	150	g

Buffered USB2.0 Monochrome 1.3MP CCD Cameras (C-series)

Mightex USB2.0 CCD cameras with frame buffers are optimized for machine-vision applications, and they can be also used for a wide variety of applications (such as industrial inspections, digital microscopy and medical imaging), which require good quality cameras that are easy to use and cost-effective. These cameras have built-in frame buffers, external trigger-in, strobe-out, and a powerful camera engine that supports multiple cameras. Monochrome cameras often exhibit 20% higher spatial resolution than their color counterparts because no pixel interpolation is necessary. Since there is no *Bayer color filter* on the sensor, monochrome cameras are more sensitive than color sensors, especially in near IR and UV regions. Frame rate can be as high as 20 fps in full resolution and up to 84 fps using ROI mode. A user-friendly GUI based application software and an SDK are provided for custom software development. In addition, a DirectShow driver and a TWAIN driver are available to easily link the cameras with users' applications. A USB command set protocol and a Linux driver are also provided for non-Windows based applications.




PERFORMANCE SPECIFICATIONS

Parameters	CGN-B013-U	CGE-B013-U	Unit
Board-level / enclosed	Board-level	Enclosed	
Number of GPIOs	4	4	
Resolution	1280x960 monochrome		
Bit	8/12		bits
CCD Chip	1/3" Sony ICX445AL, Global Shutter		
Pixel Size	3.75x3.75		μm
Active Imager Size	6.26x5.01		mm
Scanning System	Progressive		
Frame Rates (@32MHz Clock)	20 @1280 x 960, 38 @640 x 480 (2x2 Bin), 53 @424 x 320 (3x3 Bin), 66 @320 x 240 (4x4 Bin), 66 @320 x 240 (4x4 Bin2)		fps
Sub Resolutions	640 x 480 (2x2 Bin), 424 x 320 (3x3 Bin), 320 x 240 (4x4 Bin), 320 x 240 (1:4 Skip)**		
Shutter Speed (Exposure time)	0.05 ~ 200,000		ms
Hardware Gains	6 ~ 41		dB
Trigger Mode	With external trigger		
Trigger Delay	<25		μs
Strobe Out	Yes		
Lens mount	CS- mount or C-mount		
Built-in Filters	IR-cut (factory standard), IR-pass or No filter		
Power consumption	< 1.8		W
Dimension	40 x 40 x 31 (CS-mount) 40 x 40 x 36 (C-mount)	45 x 45 x 30.5 (CS-mount) 45 x 45 x 35.5 (C-mount)	mm
Weight (excluding lens)	29	115	g

Buffered USB2.0 Color 1.3MP CCD Cameras (C-Series)

Mightex USB2.0 CCD cameras with frame buffers are optimized for machine-vision applications, and they can be also used for a wide variety of applications (such as industrial inspections, digital microscopy and medical imaging) which require good quality cameras that are easy to use and cost-effective. These cameras have built-in frame buffers, external trigger-in, strobe-out, and a powerful camera engine that supports multiple cameras. Frame rate can be as high as 20 fps in full resolution and up to 84 fps using ROI mode. A user-friendly GUI based application software and an SDK are provided for custom software development. In addition, a DirectShow driver and a TWAIN driver are available to easily link the cameras with users' applications. A USB command set protocol and a Linux driver are also provided for non-Windows based applications.



PERFORMANCE SPECIFICATIONS

Parameters	CGN-C013-U	CGE-C013-U	Unit
Board-level / enclosed	Board-level	Enclosed	
Number of GPIOs	4	4	
Resolution	1280x960 Color (RGB Bayer Filter)		
Bit	8/12		bits
CCD Chip	1/3" Sony ICX445AK, Global Shutter		
Pixel Size	3.75x3.75		μm
Active Imager Size	6.26x5.01		mm
Scanning System	Progressive		
Frame Rates (@32MHz Clock)	20 @1280 x 960, 38 @640 x 480 (2x2 Bin), 53 @424 x 320 (3x3 Bin), 66 @320 x 240 (4x4 Bin), 66 @320 x 240 (4x4 Bin2)		fps
Sub Resolutions	640 x 480 (2x2 Bin), 424 x 320 (3x3 Bin), 320 x 240 (4x4 Bin), 320 x 240 (1:4 Bin2)		
Shutter Speed (Exposure time)	0.05 ~ 200,000		ms
Hardware Gains	6 ~ 41		dB
Trigger Mode	With external trigger		
Trigger Delay	<25		μs
Strobe Out	Yes		
Lens mount	CS- mount or C-mount		
Built-in Filters	IR-cut (factory standard), IR-pass or No filter		
Power consumption	< 1.8		W
Dimension	40 x 40 x 31 (CS-mount) 40 x 40 x 36 (C-mount)	45 x 45 x 30.5 (CS-mount) 45 x 45 x 35.5 (C-mount)	mm
Weight (excluding lens)	29	115	g

Buffered USB2.0 Monochrome 1.4MP CCD Cameras (C-Series)

Mightex USB2.0 cameras with frame buffers are optimized for machine-vision applications, and they can be also used for a wide variety of applications (such as industrial inspections, digital microscopy and medical imaging), which require good quality cameras that are easy to use and cost-effective. These cameras have built-in frame buffers, external trigger-in, strobe-out, and a powerful camera engine that supports multiple cameras. Monochrome cameras often exhibit 20% higher spatial resolution than their color counterparts because no pixel interpolation is necessary. Since there is no *Bayer color filter* on the sensor, monochrome cameras are more sensitive than color sensors, especially in near IR and UV regions. Frame rate can be as high as 15 fps in full resolution and up to 80 fps using ROI mode. A user-friendly GUI based application software and an SDK are provided for custom software development. In addition, a DirectShow driver and a TWAIN driver are available to easily link the cameras with users' applications. A USB command set protocol and a Linux driver are also provided for non-Windows based applications.



PERFORMANCE SPECIFICATIONS

Parameters	CCN-B013-U	CCE-B013-U	Unit
Board-level/enclosed	Board-level	Enclosed	
Number of GPIOs	4	4	
Resolution	1392x1040 monochrome		
Bit	8/12		bits
CCD Chip	1/2" Sony ICX205AL, Global shutter		
Pixel Size	4.65x4.65		μm
Active Imager Size	7.60x6.20		mm
Scanning System	Progressive		
Frame Rates (@28MHz Clock)	15 @1392 x 1040, 29 @696 x 520 (2x2 Bin), 37 @464 x 344 (3x3 Bin), 49 @348 x 256 (4x4 Bin), 49 @348 x 256 (1:4 Skip)		fps
Sub Resolutions	696 x 520(2x2 Bin), 464 x 344 (3x3 Bin), 348 x 256 (4x4 Bin), 348 x 256 (1:4 Skip)		
Shutter Speed (Exposure time)	0.05 ~ 200,000		ms
Hardware Gains	6 ~ 43		dB
Trigger Mode	With external trigger		
Trigger Delay	<25		ms
Strobe Out	Yes		
Lens mount	CS- mount or C-mount		
Built-in Filters	IR-cut (factory standard), IR-pass or No filter		
Power consumption	< 1.8		W
Dimension	89 x 64 x 34 (CS-mount) 89 x 64 x 39 (C-mount)	95 x 70 x 38.5 (CS-mount) 95 x 70 x 43.5 (C-mount)	mm
Weight (excluding lens)	80	150	g

Buffered USB2.0 Color 1.4MP CCD Cameras (C-Series)

Mightex USB 2.0 CCD cameras with frame buffers are optimized for machine-vision applications, and they can be also used for a wide variety of applications (such as industrial inspections, digital microscopy and medical imaging), which require good quality cameras that are easy to use and cost-effective. These cameras have built-in frame buffers, external trigger-in, strobe-out, and a powerful camera engine that supports multiple cameras. Frame rate can be as high as 15 fps in full resolution and up to 80 fps using ROI mode. A user-friendly GUI based application software and an SDK are provided for custom software development. In addition, a DirectShow driver and a TWAIN driver are available to easily link the cameras with users' applications. A USB command set protocol and a Linux driver are also provided for non-Windows based applications.



PERFORMANCE SPECIFICATIONS

Parameters	CCN-C013-U	CCE-C013-U	Unit
Board-level/enclosed	Board-level	Enclosed	
Number of GPIOs	4	4	
Resolution	1392x1040 Color		
Bit	8/12		bits
CCD Chip	1/2" Sony ICX205AK, Global shutter		
Pixel Size	4.65x4.65		μm
Active Imager Size	7.60x6.20		mm
Scanning System	Progressive		
Frame Rates (@28MHz Clock)	15 @1392 x 1040, 49 @348 x 256 (1:4 Skip)*, Resolution dependent @ ROI		fps
Sub Resolutions	348 x 256 (1:4 Skip)		
Shutter Speed (Exposure time)	0.05 ~ 200,000		ms
Hardware Gains	6 ~ 43		dB
Trigger Mode	With external trigger		
Trigger Delay	<25		μs
Strobe Out	Yes		
Lens mount	CS- mount or C-mount		
Built-in Filters	IR-cut (factory standard), IR-pass or No filter		
Power consumption	< 1.8		W
Dimension	89 x 64 x 34 (CS-mount) 89 x 64 x 39 (C-mount)	95 x 70 x 38.5 (CS-mount) 95 x 70 x 43.5 (C-mount)	mm
Weight (excluding lens)	80	150	g

High-Sensitivity USB2.0 Monochrome 8/12-Bit 1.4MP 2/3" CCD Camera (C-Series)

Mightex USB 2.0 CCD cameras with frame buffers are optimized for machine-vision applications, and they can be also used for a wide variety of applications (such as industrial inspections, digital microscopy and medical imaging), which require good quality cameras that are easy to use and cost-effective. These cameras have built-in frame buffers, external trigger-in, strobe-out, and a powerful camera engine that supports multiple cameras. Monochrome cameras often exhibit 20% higher spatial resolution than their color counterparts because no pixel interpolation is necessary. Without a Bayer color filter on the sensor, monochrome cameras are also more sensitive than color sensors especially in near IR and UV regions. Frame rate can be as high as 15 fps in full resolution and up to 80 fps using ROI mode. In addition, a user-friendly GUI based application software and an SDK are provided for custom software development. A USB command set protocol is also provided for non-Windows based applications.



PERFORMANCE SPECIFICATIONS

Parameters	CXN-B013-U	CXE-B013-U	Unit
Board-level/enclosed	Board-level	Enclosed	
Number of GPIOs	4	4	
Resolution	1392x1040 monochrome		
Bit	8/12		bits
CCD Chip	2/3" Sony ICX285AL, Global shutter		
Pixel Size	6.45x6.45		μm
Active Imager Size	Diagonal 11.0		mm
Scanning System	Progressive		
Frame Rates (@28MHz Clock)	15 @1392 x 1040, 29 @696 x 520 (2x2 Bin), 37 @464 x 344 (3x3 Bin), 49 @348 x 256 (4x4 Bin), 49 @348 x 256 (1:4 Skip)		fps
Sub Resolutions	696 x 520 (2x2 Bin), 464 x 344 (3x3 Bin), 348 x 256 (4x4 Bin), 348 x 256 (1:4 Skip)		
Shutter Speed (Exposure time)	0.05~200,000		ms
Hardware Gains	6~43		dB
Trigger Mode	With external trigger		
Trigger Delay	<25		us
Strobe Out	Yes		
Lens mount	CS- mount or C-mount		
Built-in Filters	IR-cut (factory standard), IR-pass or no filter		
Power consumption	< 1.8		W
Dimension	89 x 64 x 34 (CS-mount) 89 x 64 x 39 (C-mount)	95 x 70 x 38.5 (CS-mount) 95 x 70 x 43.5 (C-mount)	mm
Weight (excluding lens)	80	150	g

High-Sensitivity USB2.0 Color 8/12-Bit 1.4MP 2/3" CCD Camera (C-Series)

Mightex USB 2.0 CCD cameras with frame buffers are optimized for machine-vision applications, and they can be also used for a wide variety of applications (such as industrial inspections, digital microscopy and medical imaging), which require good quality cameras that are easy to use and cost-effective. These cameras have built-in frame buffers, external trigger-in, strobe-out, and a powerful camera engine that supports multiple cameras. Frame rate can be as high as 15 fps in full resolution and up to 80 fps using ROI mode. In addition, a user-friendly GUI based application software and an SDK are provided for custom software development. A USB command set protocol is also provided for non-Windows based applications.




PERFORMANCE SPECIFICATIONS

Parameters	CXN-C013-U	CXE-C013-U	Unit
Board-level/enclosed	Board-level	Enclosed	
Number of GPIOs	4	4	
Resolution	1392x1040 color		
Bit	8/12		bits
CCD Chip	2/3" Sony ICX285AQ, Global shutter		
Pixel Size	6.45x6.45		μm
Active Imager Size	Diagonal 11.0		mm
Scanning System	Progressive		
Frame Rates (@28MHz Clock)	15 @1392 x 1040, 49 @348 x 256 (1:4 Skip)		fps
Sub Resolutions	348 x 256 (1:4 Skip)		
Shutter Speed (Exposure time)	0.05~200,000		ms
Hardware Gains	6~43		dB
Trigger Mode	With external trigger		
Trigger Delay	<25		us
Strobe Out	Yes		
Lens mount	CS- mount or C-mount		
Built-in Filters	IR-cut (factory standard), IR-pass or No filter		
Power consumption	< 1.8		W
Dimension	89 x 64 x 34 (CS-mount) 89 x 64 x 39 (C-mount)	95 x 70 x 38.5 (CS-mount) 95 x 70 x 43.5 (C-mount)	mm
Weight (excluding lens)	80	150	g

Windowless USB2.0 Monochrome 1.3MP CMOS Camera (8-bit)

Mightex USB 2.0 monochrome cameras are designed for a wide variety of applications such as digital microscopy and medical imaging, where quality, ease of use, and cost-effectiveness are crucial. Monochrome cameras often exhibit 20% higher spatial resolution than their color counterparts because no pixel interpolation is necessary. Without a Bayer color filter on the sensor, monochrome cameras are also more sensitive than color sensors especially in near IR and UV regions. Frame rate can be as high as 22 fps in full resolution and up to 600 fps using ROI mode. In addition, a user-friendly GUI based application software and an SDK are provided for custom software development. The cameras have 4-pin GPIOs and a TWAIN driver, a DirectShow driver and a LabVIEW driver are also provided to easily link the cameras with users' applications. Window-less cameras are particularly useful for applications that are sensitive to multiple reflections due to the existence of the glass window in front of the CMOS sensor.




PERFORMANCE SPECIFICATIONS

Parameters	MCE-B013-UW and MCE-B013-UWS	Unit
Board-level/enclosed	Enclosed	
Number of GPIOs	4	
Resolution	1280x1024 monochrome	
CMOS Chip	½" (5:4) Aptina MT9M001	
Pixel Size	5.2 x 5.2	μm
Active Imager Size	6.66 x 5.32	mm
Scanning System	Progressive	
Dynamic Range	68	dB
Sensor SNR	45	dB
Gray Level	8	bit
Responsivity	2.1	V/lux-sec
Frame Rates (@48MHz Clock)	22 @1280x1024, 35 @1024x768, 50 @800x600, 70 @640x480, 200 @320x240, 300 @160x120, 450 @64x64, 600 @32x32	fps
Sub Resolutions	1024x768, 800x600, 640x480, 320x240, 160x120, 64x64, 32x32	
Shutter Speed (Exposure time)	0.04~750	ms
Hardware Gains (0x – 16x)	12	dB
Trigger Mode	With external trigger	
Lens mount	C-Mount, CS-Mount or custom-defined	
Built-in Filters	None	
Power consumption	< 1.0	W
Dimension	58 x 58 x 34 (CS-mount), 58 x 58 x 39 (C-mount)	mm
Weight (excluding lens)	150	g

Windowless USB2.0 Buffered CMOS Monochrome 1.3MP Cameras (10-bit)


Mightex USB 2.0 monochrome cameras are designed for a wide variety of applications such as digital microscopy and medical imaging, where quality, ease of use, and cost-effectiveness are crucial. Monochrome cameras often exhibit 20% higher spatial resolution than their color counterparts because no pixel interpolation is necessary. Without a Bayer color filter on the sensor, monochrome cameras are also more sensitive than color sensors especially in near IR and UV regions. Frame rate can be as high as 12 fps in full resolution and up to 150 fps using ROI mode. A user-friendly GUI based application software and an SDK are provided for custom software development. In addition, a TWAIN driver, a DirectShow driver and a LabVIEW driver are also provided to easily link the cameras with users' applications. A USB command set protocol and a Linux driver are also provided for non-windows based applications. Window-less cameras are particularly useful for applications that are sensitive to multiple reflections due to the existence of the glass window in front of the CMOS sensor.

**PERFORMANCE SPECIFICATIONS**

Parameters	BTE-B013-UW and BTE-B013-UWS	Unit
Board-level/enclosed	Enclosed	
Number of GPIOs	4	
Resolution	1280x1024 monochrome	
CMOS Chip	½" (5:4) Aptina MT9M001	
Pixel Size	5.2 x 5.2	μm
Active Imager Size	6.66 x 5.32	mm
Scanning System	Progressive	
Dynamic Range	68	dB
Sensor SNR	45	dB
Gray Level	10	bit
Responsivity	2.1	V/lux-sec
On-board Memory ("Frame Buffers")	32	MB
Frame Rates (@48MHz Clock)	12 @1280x1024, 18 @1024x768, 25 @800x600, 28 @752x480, 33 @640x480, 70 @320x240	fps
Sub Resolutions	Support arbitrary ROI, (X, Y) with X and Y multiples of 4	
Shutter Speed (Exposure time)	0.04~750	ms
Hardware Gains	0.125x ~ 8x	
Trigger Mode	With external trigger	
Lens mount	C-Mount, CS-Mount or custom-defined	
Built-in Filters	None	
Power consumption	< 1.8	W
Dimension	58 x 58 x 34 (CS-mount) 58 x 58 x 39 (C-mount)	mm
Weight (excluding lens)	150	g

3648 Pixel 16-bit CCD Line Camera with External Trigger (Board-Level)

Mightex's TCN-1304-U line camera is a cost-effective high-performance B/W board-level line camera, based on a single-line, 3648-pixel CCD chip with USB2.0 (480 Mb/s) interface. CCD line cameras have several advantages over their area-array counterparts, including high optical linear resolution that allows systems developers to use the cameras to capture two-dimensional (2-D) images by moving the object or the CCD perpendicularly to the scan line. The TCN-1304-U is a compact, board-level line-scan camera ideal for a variety of OEM applications in industry process control, optical spectroscopy and bio-medical imaging etc. Setting up the TCN-1304-U line camera is very easy, the user simply installs the latest version of the operating software onto any desktop or notebook PC and then connects the USB cable from the line camera to the PC. There is no need for installing a DAC card, or using an external power supply. A user-friendly GUI based application software, LabVIEW support and an SDK are provided. In addition a USB command set protocol and a Linux driver are available for non-Windows based applications.



"Optimized for optical spectrometers"


PERFORMANCE SPECIFICATIONS

Parameters	TCN-1304-U	Unit
Board-level/enclosed	Board-level	
CCD chip	Toshiba TCD1304DG	
Number of pixels	3648	pixels
Pixel size	8 x 200	μm
Spectral range	350 ~1000	nm
Pixel output clock	0.5	MHz
Data storage on camera	4	frames
ADC resolution	16	bits
External trigger	Yes	
Exposure time range	0.1 ~ 6,500	ms
Number of GPIOs	4 programmable I/O's	
Frame rate	138	scans/second
Host interface	USB2.0	

SDK Features	
Operation Systems	Windows 2000, XP, Vista and Windows 7
Minimum Requirement	RAM > 64M, hard disk space used > 10M
USB Port	2.0
Multiple Cameras	Supported
Device Driver	Yes
Demo Application	Yes
Library Files	Yes (DLL files and static library file)
Example Codes	Yes (VC++ and Delphi)
Frame Attributes	Exposure Time, Time Stamp, Trigger Event Count, Over-Exposure Detection.

3648 Pixel 16-bit CCD Line Camera with External Trigger (Enclosed)

Mightex's TCE-1304-U line camera is a cost-effective high-performance B/W enclosed line camera, based on a single-line, 3648-pixel CCD chip with USB2.0 (480 Mb/s) interface. CCD line cameras have several advantages over their area-array counterparts, including high optical linear resolution that allows systems developers to use the cameras to capture two-dimensional (2-D) images by moving the object or the CCD perpendicularly to the scan line. TCE-1304-U is a compact line-scan camera ideal for a variety of OEM applications such as industry process control, optical spectroscopy and bio-medical imaging etc. A **windowless** version TCE-1304-UW is available which has a broader spectral response of 200 ~ 1000. Setting up the TCE-1304-U line camera is very easy: one simply needs to install the camera's application software into any PC, and then connect the line camera to the PC using a USB cable. There is no need to install a DAC card or to use an external power supply. A user-friendly GUI based application software, LabVIEW support and an SDK are provided. In addition a USB command set protocol and a Linux driver are available for non-Windows based applications.



PERFORMANCE SPECIFICATIONS

Parameters	TCE-1304-U and TCE-1304-UW	Unit
Board-level/enclosed	Enclosed	
CCD chip	Toshiba TCD1304DG	
Number of pixels	3648	pixels
Pixel size	8 x 200	μm
Spectral range	TCE-1304-U: 350 ~1000 TCE-1304-UW: 200 ~ 1000	nm nm
Pixel output clock	0.5	MHz
Data storage on camera	4	frames
ADC resolution	16	bits
External trigger	Yes	
Exposure time range	0.1 ~ 6,500	ms
Number of GPIOs	4 programmable I/O's	
Frame rate	138	scans/second
Dimensions (Length x Width x Height)	78 x 42 x 57	mm
Weight	200	g

SDK Features	
Operation systems	Windows 2000, XP, Vista and Windows 7
Minimum requirement	RAM > 64M, Hard Disk Space Used > 10M
USB port	2.0
Multiple cameras	Supported
Device driver	Yes
Demo application	Yes
Library files	Yes (DLL files and static library file)
Example codes	Yes (VC++, CSharp, Visual Basic and Delphi)
Frame Attributes*	Exposure Time, Time Stamp, Trigger Event Count, Over-Exposure Detection.

High-Speed 2048-Pixel 12-bit CCD Line Camera with External Trigger

Mightex's TCN/TCE-1209-U line camera is a cost-effective high-performance monochrome line CCD camera, based on a single-line, 2048-pixel CCD chip with USB2.0 (480 Mb/s) interface. CCD line cameras have several advantages over their area-array counterparts, including high optical linear resolution that allows systems developers to use the cameras to capture two-dimensional (2-D) images by moving the object or the CCD perpendicularly to the scan line. The TCN/TCE-1209-U is a compact linear camera ideal for a variety of OEM applications in industry process control, optics, biology, spectroscopy and reflection imaging. Setting up the TCN/TCE-1209-U line camera is very easy: the user simply installs the application software (included in the shipping package) and then connects the line camera to the PC using a USB cable. There is no need to install a DAC card, and no external power supply is required.



PERFORMANCE SPECIFICATIONS

Parameter	TCN-1209-U or TCE-1209-U
CCD	Toshiba TCD1209
Number of Pixels	2048
Pixel Size	14 μ m x 14 μ m
Spectral range	350~1000nm
Pixel Output Clock	8 MHz
Frame buffers on camera	64 frames
ADC resolution	12 bits
External Trigger	Yes
Exposure Time Range	0.3 ~ 3,277ms
GPIO	Yes (4 programmable I/Os)
Frame Rate	3,300 scans/second*
Lens mount (optional)	F mount
Host Interface	USB 2.0

SDK Features	
Operation Systems	Windows 2000, XP, Vista, Windows 7
Minimum Requirement	RAM > 64M; Available Hard Disk Space > 10M
USB Port	USB2.0
Multiple Cameras	Supported
Plug & Play	Not recommended
Device Driver	Yes
Demo Application	Yes
Library Files	Yes (DLL files and Static Library file)
Example Codes	Yes (VC++ and Delphi)
Frame Attributes	Exposure Time, Time Stamp, Trigger Event Count, Over Exposure Detect.

(This page is left blank intentionally)

Table of Content

Manual Controlled

Manual & Analog-Input Controlled Universal LED Controllers with Current Display.....30

Manual and Analog-Input Controlled Universal 2-Channel LED Controllers31

Computer Controlled

Compact Universal 1- and 2- Channel LED Controllers32

Universal 2- and 4-Channel LED Controllers with External Triggers33

High-Precision Universal 2- and 4-Channel LED Controllers with External Triggers34

High-Current Universal 2-Channel LED Controllers with External Triggers.....35

Universal 12- and 16- Channel LED Controllers.....36

Manual & Computer (“Dual Mode”) Controlled

High-Current Dual-Mode (Manual/Software) LED Controller.....37

Manual & Analog-Input Controlled Universal LED Controller with Current Display

Mightex's SLB-1200-1 universal LED driver is designed for driving a broad range of LED light sources. The LED driver has two operational modes:

- 1) **Manual Knob Control Mode:** the output current can be adjusted manually;
- 2) **Analog Input Control Mode:** the output current can be controlled via a 0~5V analog input signal.

The control mode is set via a DIP switch, and the factory default setting is "Manual Knob Control Mode". The driver also has a Maximum Current Setting DIP Switch, which allows user to set the maximum current to 350mA, 750mA or 1200mA, whichever applicable. The factory default setting is 350mA. When the Maximum Current Setting DIP Switch is set at a smaller value (e.g. 350mA), the LED driver has a finer resolution for the output current. When the driver is set to "Analog Input Control Mode", the output current is proportional to the voltage of the analog input signal.



ELECTRICAL SPECIFICATION

Parameters	SLA-1200-1	Unit
Number of Channels	1	
Power Supply Input Voltage (V_{dc})	9 ~ 24	V
Maximum Output Voltage (V_{max})	$V_{dc} - 3.0$	V
Maximum Per Channel Output Current (I_{max})*	1200	mA
Maximum Per Channel Output Power (P_{max})**	10	W

* The maximum output current can be set to 1200, 750 or 350 mA using the DIP switch.

** If the channel output voltage is V_d and the output current is I_d , they must simultaneously satisfy the following conditions:

(1) $V_d \leq V_{max}$; (2) $I_d \leq I_{max}$; and (3) $V_d * I_d \leq P_{max}$.

CHANNEL I/O PIN DEFINITION

Each channel has four pins, defined as follows:

Pin Label	LED+	LED-	Analog Signal	Analog Input GND
Description	LED Anode	LED Cathode	0~5V Analog Input	Analog Input Ground

Manual and Analog-Input Controlled Universal 2-Channel LED Controllers

Mightex's SLA-series 2-channel universal LED controllers are designed to drive a broad range of LED light sources. These LED drivers have two operational modes:

- 1) **Manual Knob Control Mode:** the current output of each channel can be adjusted manually;
- 2) **Analog Input Control Mode:** the current output of each channel can be controlled via 0 ~ 5V analog input



The control mode is set via a DIP switch, and the factory default setting is "Manual Knob Control Mode". The drivers also have a Maximum Current Setting DIP Switch, which allows user to set the maximum current to 350mA, 750mA or 1200mA, whichever applicable. The factory default setting is 350mA. When the Maximum Current Setting DIP Switch is set at a smaller value (e.g. 350mA), the LED driver has a finer resolution for the output current. When the driver is set to "Analog Input Control Mode", the output current is proportional to the voltage of the analog input signal. For the 2-channel models, the operational mode and the current limit of each channel can be set independently from each other.

ELECTRICAL SPECIFICATION

Parameters	SLA-0350-2	SLA-0750-2	SLA-1200-2	Unit
Number of Channels	2	2	2	
Power Supply Input Voltage (V_{dc})	9 ~ 24	9 ~ 24	9 ~ 24	V
Maximum Output Voltage (V_{max}) ¹	<21	<21	<21	V
Maximum Per Channel Output Current (I_{max}) ²	350	750	1200	mA
Maximum Per Channel Output Power (P_{max}) ³	10	10	10	W

1. Maximum Output Voltage is 3V less than the Power Supply Input Voltage, i.e. $V_{max} = V_{dc} - 3V$. For instance, with a Power Supply Input Voltage of $V_{dc}=24V$, the Maximum Output Voltage V_{max} would be $(V_{dc}-3V)=21V$; and

2. If the channel output voltage is V_d and the output current is I_d , they must simultaneously satisfy: (1) $V_d \leq V_{max}$; (2) $I_d \leq I_{max}$; and (3) $V_d * I_d \leq P_{max}$.

3. If the channel output voltage is V_d and the output current is I_d , they must simultaneously satisfy: (1) $V_d \leq V_{max}$; (2) $I_d \leq I_{max}$; and (3) $V_d * I_d \leq P_{max}$.

CHANNEL I/O PIN DEFINITION

Each channel has four pins, defined as follows:

Pin Label	LED+	LED-	Analog Signal	Analog Input GND
Description	LED Anode	LED Cathode	0~5V Analog Input	Analog Input Ground

Compact Universal 1- and 2- Channel LED Controllers

Mightex's Compact 1- and 2- Channel Computer-Controlled Universal LED Drivers are designed to drive a broad range of LED light sources. Each unit comes with a powerful PC-based software with a user-friendly GUI, which enables users to drive LEDs without the need to write any code. In addition, a full-featured SDK is provided, in order for users to write their own software and to integrate Mightex's LED drivers into their own systems. Furthermore, the drivers have a built-in protection feature, allowing users to limit LED driving current and voltage.

Each channel can be individually configured to work under one of the following two modes:

- 1) **Normal Mode (or DC Mode):** The output current is a constant, which can be adjusted (using software) from 0 mA to 1,000 mA, through the USB interface;
- 2) **Strobe Mode:** A Pulse-Width-Modulated (or PWM) periodic strobe pattern is output from the channel, which can be turned on by a software trigger. The strobe pattern may last indefinitely or for a preset number of cycles. The frequency of the PWM strobe can be up to 500Hz. In addition, each channel can be individually **DISABLED** and **ENABLED**. No voltage or current is output from a DISABLED channel.



ELECTRICAL SPECIFICATION

Parameters	SLC-MA01-U	SLC-MA02-U	Unit
Number of Channels	1	2	
Power Supply Input Voltage (V_{dc})	9 ~ 24		V
Maximum Output Voltage (V_{max}) ¹	<21		V
Maximum Per Channel Output Current (I_{max})	1,000		mA
Maximum Per Channel Output Power (P_{max}) ²	10		W
Output Current Resolution	1		mA
Output Current Accuracy	± 5 mA or $\pm 1.0\%$, whichever is larger		mA
Output Current Repeatability	± 2 mA or $\pm 0.5\%$, whichever is larger		mA
PWM Timing Resolution ³	100		ms
PWM Timing Minimum Step Size ³	1,000		ms
Interface	USB		

1. Maximum Output Voltage is 3V less than the Power Supply Input Voltage, i.e. $V_{max} = V_{dc} - 3V$. For instance, with a Power Supply Input Voltage of $V_{dc}=24V$, the Maximum Output Voltage V_{max} would be $(V_{dc}-3V)=21V$;

2. If the channel output voltage is V_d and the output current is I_d , they must simultaneously satisfy: (1) $V_d \leq V_{max}$; (2) $I_d \leq I_{max}$; and (3) $V_d * I_d \leq P_{max}$; and

3. Each period of a PWM square wave comprises of ON time and OFF time, i.e. two (2) 'steps'. The minimum value for each step is 1000ms, and the minimum increment is 100ms.

CHANNEL I/O PIN DEFINITION

Each Channel has two pins, defined as following:

Label	LED+	LED-
Description	LED Anode	LED Cathode

Universal 2- and 4-Channel LED Controllers with External Triggers

Mightex has developed a series of computer-controllable, multi-channel, universal LED controllers, which can be used to drive any type of LED in any of the three (3) modes: 'NORMAL' (or 'constant current'), 'STROBE', and/or external 'TRIGGER' mode. Each unit comes with PC-based software with a user-friendly GUI, which enables users to drive LEDs without the need to write any code. In addition, a powerful SDK is provided, in order for users to write their own software and to integrate Mightex's LED drivers into their own systems. Furthermore, the drivers have a built-in security feature, allowing users to limit LED driving current and voltage.



ELECTRICAL SPECIFICATION

Parameters	SLC-AA02-US SLC-AA04-US	SLC-AV02-US SLC-AV04-US	SLC-SA02-US SLC-SA04-US	SLC-SV02-US SLC-SV04-US	Unit
Power Supply Input Voltage, V(dc)	9 ~ 24				V
Power Supply Input Current	< 4000				mA
Per Channel Driving Voltage (max)	11.5V (with a 12V DC input) or 23.5V (with a 24V DC input)				V
Per Channel Driving Current	0 ~ 1000 ("NORMAL" MODE)				mA
	0 ~ 3500 ("STROBE" or "TRIGGER" MODE)				mA
Output Current Resolution	1				mA
Output Current Linearity	+/-4 (or +/-0.5%)				mA
Output Current Repeatability	+/-1 (or +/-0.2%)				mA
Trigger Input High Level	4.5 ~ 10.0				V
Trigger Input Low Level	0.8 (Max.)				V
Forward Voltage Monitoring Capability	No	Yes	No	Yes	-
Forward Voltage Monitoring Accuracy	N.A.	+/-10	N.A.	+/-10	mV

TIMING SPECIFICATION:

Parameters	SLC-AA02-US SLC-AA04-US	SLC-AV02-US SLC-AV04-US	SLC-SA02-US SLC-SA04-US	SLC-SV02-US SLC-SV04-US	Unit
Timing Resolution	20				μs
# of Data Points for Waveform Definition	128		2		-
Trigger Pulse Width	100 (Minimum)				μs
Max Trigger Delay	25				μs

High-Precision Universal 2- and 4-Channel LED Controllers with External Triggers

Mightex has developed a series of computer-controllable, multi-channel, universal LED controllers, which can be used to drive any type of LED in any of the three (3) modes: 'NORMAL' (or 'constant current'), 'STROBE', and/or external 'TRIGGER' mode. Each unit comes with PC-based software with a user-friendly GUI, which enables users to drive LEDs without the need to write any code. In addition, a powerful SDK is provided, in order for users to write their own software and to integrate Mightex's LED drivers into their own systems. Furthermore, the drivers have a built-in security feature, allowing users to limit LED driving current and voltage.



ELECTRICAL SPECIFICATION

Parameters	SLC-FA02-US SLC-FA04-US	SLC-FV02-US SLC-FV04-US	SLC-XA02-US SLC-XA04-US	SLC-XV02-US SLC-XV04-US	Unit
Power Supply Input Voltage, V(dc)	9 ~ 24				V
Power Supply Input Current	< 4000				mA
Per Channel Driving Voltage (max)	11.5V (with a 12V DC input) or 23.5V (with a 24V DC input)				V
Per Channel Driving Current	0 ~ 100 ("NORMAL" MODE)				mA
	0 ~ 350 ("STROBE" or "TRIGGER" MODE)				mA
Output Current Resolution	0.1				mA
Output Current Linearity	+/-4 (or +/-0.5%)				mA
Output Current Repeatability	+/-1 (or +/-0.2%)				mA
Trigger Input High Level	4.5 ~ 10.0				V
Trigger Input Low Level	0.8 (Max.)				V
Forward Voltage Monitoring Capability	No	Yes	No	Yes	-
Forward Voltage Monitoring Accuracy	N.A.	+/-10	N.A.	+/-10	mV

TIMING SPECIFICATION:

Parameters	SLC-FA02-US SLC-FA04-US	SLC-FV02-US SLC-FV04-US	SLC-XA02-US SLC-XA04-US	SLC-XV02-US SLC-XV04-US	Unit
Timing Resolution	20				μs
# of Data Points for Waveform Definition	2		128		-
Trigger Pulse Width	100 (Minimum)				μs
Max Trigger Delay	25				μs

High-Current Universal 2-Channel LED Controllers with External Triggers

Mightex has developed a series of computer-controllable, multi-channel, universal LED controllers, which can be used to drive any type of LED in any of the three (3) modes: 'NORMAL' (or 'constant current'), 'STROBE', and/or external 'TRIGGER' mode. Each unit comes with PC-based software with a user-friendly GUI, which enables users to drive LEDs without the need to write any code. In addition, a powerful SDK is provided, in order for users to write their own software and to integrate Mightex's LED drivers into their own systems. Furthermore, the drivers have a built-in security feature, allowing users to limit LED driving current and voltage.



ELECTRICAL SPECIFICATION

Parameters	SLC-HA02-US	SLC-HV02-US	Unit
Power Supply Input Voltage, V(dc)	9 ~ 12		V
Power Supply Input Current	4,000		mA
Per Channel Driving Voltage (max)	Up to 11.5		V
Maximum Per Channel Driving Current	2,000 ("NORMAL" MODE)		mA
	3,500 ("STROBE" or "TRIGGER" MODE)		mA
Output Current Resolution	1		mA
Output Current Linearity	+/-4 (or +/-0.5%)		mA
Output Current Repeatability	+/-1 (or +/-0.2%)		mA
Trigger Input High Level	3.3 ~ 10.0		V
Trigger Input Low Level	0.8 (Max.)		V
Forward Voltage Monitoring Accuracy	N/A	+/-10	mV

TIMING SPECIFICATION:

Parameters	SLC-HA02-US	SLC-HV02-US	Unit
Timing Resolution	20		μs
# of Data Points for Waveform Defini-	2		
Trigger Pulse Width	100 (Minimum)		μs
Max Trigger Delay	25		μs

Universal 12- and 16- Channel LED Controllers

Mightex's 12- and 16- Channel Computer-Controlled Universal LED Drivers are designed to drive a broad range of LED light sources. Each unit comes with a powerful PC-based software with a user-friendly GUI, which enables users to drive LEDs without the need to write any code. In addition, a full-featured SDK is provided, in order for users to write their own software and to integrate Mightex's LED drivers into their own systems. Furthermore, the drivers have a built-in protection feature, allowing users to limit LED driving current and voltage.



Each channel can be individually configured to work under one of the following two modes:

Normal Mode (or DC Mode): The output current is a constant, which can be adjusted (using software) from 0 mA to 1,000 mA, through the USB interface;

Strobe Mode: A Pulse-Width-Modulated (or PWM) periodic strobe pattern is output from the channel, which can be turned on by a software trigger. The strobe pattern may last indefinitely or for a preset number of cycles. The frequency of

In addition, each channel can be individually **DISABLED** and **ENABLED**. No voltage or current is output from a DISABLED channel.

ELECTRICAL SPECIFICATION

Parameters	SLC-CA12-U SLC-CA12-S	SLC-CA16-U SLC-CA16-S	SLC-MA12-U SLC-MA12-S	SLC-MA16-U SLC-MA16-S	Unit
Number of Channels	12	16	12	16	
Power Supply Input Voltage (V_{dc})	9 ~ 24		9 ~ 24		V
Maximum Output Voltage (V_{max}) ¹	<21		<21		V
Maximum Per Channel Output Current (I_{max})	1,000		1,000		mA
Maximum Per Channel Output Power (P_{max}) ²	10		10		W
Output Current Resolution	5		1		mA
Output Current Accuracy	±10 mA or ±2.0%, whichever is larger		±5 mA or ±1.0%, whichever is larger		mA
Output Current Repeatability	±5 mA or ±1%, whichever is larger		±2 mA or ±0.5%, whichever is larger		mA
PWM Timing Resolution ³	100		100		μs
PWM Timing Minimum Step Size ³	1,000		1,000		μs
Interface	USB (-U) or RS232 (-S)				

1. Maximum Output Voltage is 3V less than the Power Supply Input Voltage, i.e. $V_{max} = V_{dc} - 3V$. For instance, with a Power Supply Input Voltage of $V_{dc}=24V$, the Maximum Output Voltage V_{max} would be $(V_{dc}-3V)=21V$;

2. If the channel output voltage is V_d and the output current is I_d , they must simultaneously satisfy: (1) $V_d \leq V_{max}$; (2) $I_d \leq I_{max}$; and (3) $V_d * I_d \leq P_{max}$; and

3. Each period of a PWM square wave comprises of ON time and OFF time, i.e. two (2) 'steps'. The minimum value for each step is 1000ms, and the minimum increment is 100ms.

CHANNEL I/O PIN DEFINITION

Each channel has two pins, defined as following:

Label	LED+	LED-
Description	LED Anode	LED Cathode

High-Current Dual-Mode (Manual/Software) LED Controller

Mightex's SLB-HA01-U LED controller is a constant-current source, and it is designed as a universal LED controller capable of driving a broad range of LED's. The device has two control modes:

- **Software Control Mode:** The device can be controlled via USB interface by software on a host (e.g. a PC).
- **Manual Control Mode:** When the device is set to NORMAL mode, the device can be manually controlled via a knob and three buttons on its front panel.

The SLB-HA01-U LED controller has one (1) output channel. It has a USB interface through which the device can be controlled by a host (e.g. a PC). In addition, it also has a turning knob, three buttons, and a 2x16 LCD display on its front panel, through which user can control the output current manually. In addition, one can set the maximum output current (I_{max}), manually or by software, and then adjust the actual output current ranging from 0 to I_{max} via the knob or via software.



ELECTRICAL SPECIFICATION

Parameters	SLB-HA01-U	Unit
Power Supply Input Voltage,	9 ~ 24	V(dc)
Max Output Voltage	V(dc) - 1.0	V
Channel Driving Current	0~2000 (in CW mode)	mA
	0~3500 (in Strobe/Trigger mode)	mA
Current Resolution	1 (in 0~3500mA range)	mA
Current Accuracy	+/-4 or +/-0.5% (whichever is larger)	mA
Current Repeatability	+/-1 or +/-0.2% (whichever is larger)	mA
Trigger Input High Level	4.5~10.0	V
Trigger Input Low Level 0~0.8 V	0~0.8	V

CHANNEL I/O PIN DEFINITION:

The channel has four pins as follows:

Pin Label	LED +	LED -	Trigger +	Trigger -
Description	LED Anode	LED Cathode	External Trigger +	External Trigger -

(This page is left blank intentionally)

Table of Content

Fiber-Coupled LED's	40
Light-Guide-Coupled LED's	41
LED Collimator Sources	42
Low-cost LED Spot Lights	43
Precision LED Spot Lights	45
Spectrum Synthesizing Sources ("Cubic S")	46

Industrial
Cameras

LED
Controllers

Light
Sources

Microscope
Tools

Spectrometers

High Power Fiber-Coupled LED Light Sources (UV, VIS and NIR)

Mightex FCS-series fiber-coupled LED light sources employ the latest high-power LED technologies and a proprietary coupling optics to achieve maximum optical output power. Optical output is coupled into a fiber through a standard SMA fiber adaptor port (SMA fiber patch cords are sold separately). FCS series also features a locking electrical connector for secured connection. FCS series are designed as a universal light source for general lab use and OEM applications. The one-piece machined housing features multiple mounting holes. All Mightex LED drivers such as the SLC series or other LED drivers and current sources can be used to drive the light sources.



PERFORMANCE SPECIFICATIONS

Wavelength Code	Description	Peak Wavelength (nm)	I _{op} (mA)	V _{op} (V)	Radiant Flux ¹ (mW)
0000	Cool White	5500K	1000	3.9	3.2
0001	Warm White	3400K	1000	3.9	3.2
0365	UV 365nm	365	500	3.8	3.1
0385	UV 385nm	385	500	3.8	3.4
0400	UV 400nm	400	350	3.5	0.8
0455	Royal Blue	455	1000	3.9	5.6
0470	Blue	470	1000	3.9	4.3
0505	Cyan	505	1000	3.9	1.0
0530	Green	530	1000	3.9	2.0
0590	Amber	590	1000	3.2	1.2
0617	Red-Orange	617	1000	3.0	6.5
0625	Red	625	1000	3.0	6.5
0657	Deep Red	657	350	2.4	1.0
0740	NIR	740	350	2.1	1.5
0850	NIR	850	1000	2.1	4.0
0940	NIR	940	1000	2.4	4.0

Note: * - With 400μm-core 0.22 numerical aperture (NA) fiber. Output optical power scales approximately linearly with fiber core area and NA².

High Power Light-Guide Coupled LED Light Sources (UV, VIS and NIR)

Mightex GCS-series high-power LED sources are designed for high-efficiency coupling of LED light into a liquid light guide (LLG) or a fiber optic bundle. Virtually all light guides with core diameter ranging from 3mm to 8mm can be used with the GCS-series light source. Please note that light guides and adapters are sold separately. GCS series also features a locking electrical connector for secured connection. GCS series are designed as a universal light source for general lab use and OEM applications. All Mightex's LED drivers such as the SLC series or other LED drivers and current sources can be used to drive the GCS-series light sources. The one-piece machined aluminum alloy housing features integrated heatsinks and multiple mounting holes.



Note: Liquid light guide not included.

PERFORMANCE SPECIFICATIONS

Wavelength Code	Description	Peak Wavelength (nm)	I _{op} (mA)	V _{op} (V)	Radiant Flux ¹ (mW)
0000	Cool White	5500K	1000	3.9	95
0001	Warm White	3400K	1000	3.9	95
0365	UV 365nm	365	500	3.8	90
0385	UV 385nm	385	500	3.8	100
0455	Royal Blue	455	1000	3.9	165
0470	Blue	470	1000	3.9	130
0505	Cyan	505	1000	3.9	30
0530	Green	530	1000	3.9	60
0590	Amber	590	1000	3.2	35
0617	Red-Orange	617	1000	3.0	200
0625	Red	625	1000	3.0	200
0850	NIR	850	1000	2.1	125
0940	NIR	940	1000	2.4	125

Note: 1 - Measured at exiting end of a 3mm-core 0.59 numerical aperture (NA) liquid light guide.

High Power LED Collimator Sources

A LED collimator consists of a collimating lens and a LED emitter. The LED emitter is placed at the focal plane of the collimating lens which images the LED emitter into infinity. Mightex LED collimators employ a high-NA aspherical collimating lens for precision collimation and high light throughput.

LED collimators can be used as the light source in an illumination system. For example LED collimators can replace the standard lamp assembly in a microscope to provide stable, intense, and fast-modulated illumination. Typically the rest of the illumination optics will image the LED emitter onto the pupil of the imaging optics where as the lens aperture on the collimator where intensity is uniform will be projected onto the object. In other illumination applications similar arrangement should be made to produce uniform and efficient illumination.

The LED emitters are mounted directly onto the metal base of the collimator which also features an integrated heat sink. This configuration minimizes thermal resistance between the LED emitter and the heat sink resulting in better heat dissipation. The collimating lens can be adjusted if needed for precise collimation. A locking ring fixes the lens position after adjustment. The collimators have been pre-adjusted in the factory.

The LED collimators include a 1-meter cable with two bare-wire terminals at the end. The light sources can be driven by Mightex LED controllers or other LED controllers and current sources. An optional focusing module can be mounted on the front of the LED collimator to focus light into a tight spot which is an image of the LED emitter. One of the applications with the focusing module is coupling LED light into a fiber or a light guide.



PERFORMANCE SPECIFICATIONS

Part Number	Wavelength	Diverging Angle (deg.)	I _{op} (mA)	V _{op} (V)	Output Power (mW)
LCS-0365-02-22	365	+/-1.7	500	3.8	80
LCS-0385-02-22	385	+/-1.7	500	3.8	100
LCS-0400-01-22	400	+/-2.5	350	3.5	100
LCS-0455-03-22	455	+/-1.7	1000	3.9	280
LCS-0455-05-22	455	+/-5.0	700	6.8	350
LCS-0470-03-22	470	+/-1.7	1000	3.9	200
LCS-0470-05-22	470	+/-5.0	700	6.8	250
LCS-0505-03-22	505	+/-1.7	1000	3.9	135
LCS-0505-05-22	505	+/-5.0	700	6.8	170
LCS-0530-03-22	530	+/-1.7	1000	3.9	100
LCS-0530-05-22	530	+/-5.0	700	6.8	125
LCS-0590-03-22	590	+/-1.7	1000	3.9	65
LCS-0617-03-22	617	+/-1.7	1000	3.9	280
LCS-0625-03-22	625	+/-1.7	1000	3.9	280
LCS-0657-01-22	657	+/-2.5	350	2.4	100
LCS-0740-01-22	740	+/-2.5	350	2.2	65
LCS-0850-02-22	850	+/-2.5	1000	2.1	175
LCS-0940-01-22	940	+/-1.7	700	1.5	100
LCS-5500-03-22	cool white 5500K	+/-1.7	1000	3.9	170
LCS-4000-03-22	warm white 4000K	+/-1.7	1000	3.9	180

Low-Cost LED Spot Lights

High-power LEDs are a new class of light sources that have numerous applications in industry, consumer, medical, and scientific fields. However, engineers and scientists often find that they have to design and fabricate heat sinks and optical mounts, and solder electrical contacts before they can light up a high-power LED. Mightex provides ready-to-use high-power LED light sources with integrated heat sinks and mounted collecting optics. Sirius™ compact high-power light sources are designed as a universal light source for general lab use and OEM applications. The mechanical housing features multiple mounting holes compatible to common opto-mechanical mounts. Sirius™ light sources can be driven by Mightex's SLC-series multi-channel LED drivers or other LED drivers and



PERFORMANCE SPECIFICATIONS

Table 1. LED Emitter Specifications



LED Code	Description	Peak Wavelength (nm)	I _{op} (mA)	V _{op} (V)	Luminous Flux (lm)
0100	1W White LED	5500K	350	3.4	45
0109	1W UV	395-410	350	3.5	180 mW
0101	1W Royal Blue	455	350	3.4	220 mW
0102	1W Blue	470	350	3.4	16
0103	1W Cyan	505	350	3.4	45
0104	1W Green	530	350	3.4	53
0105	1W Amber	590	350	3.0	42
0106	1W Red-Orange	617	350	3.0	55
0107	1W Red	625	350	3.0	44
0110	1W Deep Red	657	350	2.4	170mW
0111	1W 740nm	740	350	2.2	120mW
0112	1W 940nm	940	700	1.5	180mW
0300	3W White LED	5500K	1000	3.9	80
0309	3W UV	395-410	700	3.5	350 mW
0301	3W Royal Blue	455	1000	3.9	450 mW
0302	3W Blue	470	1000	3.9	30
0303	3W Cyan	505	1000	3.9	80
0304	3W Green	530	1000	3.9	80
0305	3W Amber	590	1400	3.0	110
0306	3W Red-Orange	617	1400	3.0	190
0307	3W Red	625	1400	3.0	140
0208	2.4W Near Infrared	850	1000	1.8~2.4	375 mW
0501	5W Royal Blue	455	1000	6.8	700 mW
0502	5W Blue	470	700	6.8	48
0503	5W Cyan	505	1000	6.8	160
0504	5W Green	530	1000	6.8	160

Low-cost LED Spot Lights (Cont'd)

Table 2. Collecting Optics Specifications

LENS Code	Description	Clear Aperture (mm)	X- Half Angle (degree)	Y- Half Angle (degree)	Efficiency (%)
A	Narrow Beam	19	5	5	85
B	Medium Beam	19	15	15	85
C	Wide Beam	19	25	25	85
D	Oval Beam	19	5	20	85
E ¹	Fiber Bundle Coupling	19	7mm full field	7mm full field	85

Part Number and Ordering Information

SLS -  - 

LED Code Lens Code

For example, SLS-0104-B is a light source with a 1W green emitter and 15-degree collecting optics.

With a world-class OEM design team, Mightex offers a broad range of customized solutions in order to meet individual customer's unique requirements. Please call 1-416-840 4991 or email


Precision LED Spot Lights

Mightex precision LED spotlight consists of a state-of-the-art high-power LED emitter and a proprietary high-NA multi-element aspherical optical system. The result is a high-power, uniform illumination spot with a highly-delineated edge.

Mightex PLS- series precision LED spotlights are a general-purpose light source that can be used where uniform and high-intensity illumination is required. The projection lens at the front of the spotlight can be slid and locked to focus the illumination pattern at different working distances. With the standard projection lens spot diameter is linearly proportional to working distance.

The LED emitters are mounted directly on the metal base of the light source which also features an integrated heatsink. This configuration minimizes thermal resistance between the LED emitter and the heatsink resulting in better heat dissipation.

The precision LED spotlight includes a 2-meter cable with two bare-wire terminals at the end. The light sources can be driven by Mightex's LED controllers, or other LED controllers and current sources.



PERFORMANCE SPECIFICATIONS

Part Number	Wavelength (nm)	I _{op} (mA)	V _{op} (V)	Output Power (mW)
PLS-0365-030-S	365	500	3.8	50
PLS-0385-030-S	385	500	3.8	50
PLS-0400-030-S	400	350	3.5	50
PLS-0455-030-S	455	1000	3.9	150
PLS-0470-030-S	470	1000	3.9	110
PLS-0505-030-S	505	1000	3.9	65
PLS-0530-030-S	530	1000	3.9	50
PLS-0590-030-S	590	1000	3.9	35
PLS-0617-030-S	617	1000	3.9	150
PLS-0625-030-S	625	1000	3.9	150
PLS-0657-030-S	657	350	2.4	50
PLS-0740-030-S	740	350	2.2	35
PLS-0850-030-S	850	1000	2.1	85
PLS-0940-030-S	940	700	1.5	50
PLS-5500-030-S	cool white 5500K	1000	3.9	85
PLS-4000-030-S	warm white 4000K	1000	3.9	85

Spectrum Synthesizing Sources (“Cubic-S”)

Many applications, such as microscopy, multi-spectral imaging, LCD display characterization, detector calibration, and color measurement etc., can benefit tremendously from a light source with an arbitrarily programmable optical spectrum. As apposed to conventional tunable lasers, which only produce a single wavelength at a time, Mightex Optical Spectrum Synthesizing Source (“Cubic-S”) is capable of dynamically generating any custom-specified target spectrum across the entire spectral range.



PRODUCT DESCRIPTION

As illustrated in Figure-1, Mightex Cubic-S is composed of the following modules: (1) An Optical Core Module; (2) A Controller Electronic Module; (3) An optional Spectrum Monitor; and (4) Control Software in the computer.




Figure 1 - Schematic diagram of Cubic-S.

Mightex's Cubic-S is controlled through a powerful PC-based application software with GUI. A software development kit (SDK) is also included, making it easy to integrate the Cubic-S into customers' own applications.

For more details regarding this product please call 1-416-840-4991 or email sales@mightex.com.

Table of Content

MLS Light Sources for Microscopes	48
Low-Cost LED Based Epi-Fluorescence Microscopy Solution.....	49
FluoFlux Fluorescence Illuminators for stereo-Microscopes	50
FilterReader Spectrophotometer	51

MLS Light Sources for Microscopes

Compatible Microscopes:

- *Leica DMI,*
- *Nikon Eclipse Microscopes,*
- *Olympus IX & BX Microscopes,*
- *Zeiss Axioskop*

Mightex microscopy light sources are collimated LED light sources designed for research microscopes. These light sources are ready to be mounted on existing lamp ports of various microscopes to replace traditional light sources such as incandescent lamps and arc lamps. MLS-series of microscopy light sources employs an aspherical collimating optics to generate uniform illuminations with high intensity at the microscope sample plane.

LEDs have many desirable features for microscopy applications. For example LED intensity can be directly modulated through driving current. LEDs also have fast response time, in some cases as short as 10ns. The output spectrum of a white LED changes little with driving current eliminating the needs for neutral density filters. Single-color LEDs have high spectral density making them suitable for applications such as multi-spectrum microscopy. With the state-of-the-art high-brightness technology spectral density of LEDs have reached or exceeded that of some short-arc lamps. LEDs may also be pulsed at current much higher than normal rating, resulting in higher output. Another feature of LEDs is their high stability and repeatability of output intensity compared to that of short-arc lamps.




PERFORMANCE SPECIFICATIONS

Part Number	Wavelength (nm)	Maximum Driving Current (mA)	Forward Voltage (V)	Output Power (mW)
MLS-0365-XXX	365	500	3.8	80
MLS-0385-XXX	385	500	3.8	100
MLS-0400-XXX	400	350	3.5	100
MLS-0455-XXX	455	1000	3.9	280
MLS-0470-XXX	470	1000	3.9	200
MLS-0505-XXX	505	1000	3.9	135
MLS-0530-XXX	530	1000	3.9	100
MLS-0590-XXX	590	1000	3.9	65
MLS-0617-XXX	617	1000	3.9	280
MLS-0625-XXX	625	1000	3.9	280
MLS-0657-XXX	657	350	2.4	100
MLS-0740-XXX	740	350	2.2	65
MLS-0850-XXX	850	1000	2.1	175
MLS-0940-XXX	940	700	1.5	100
MLS-5500-XXX	cool white 5500K	1000	3.9	170
MLS-4000-XXX	warm white 4000K	1000	3.9	180

Low-Cost LED Based Epi-Fluorescence Microscopy Solution

Episcopic(EPI) or reflective/incident illumination is commonly used in biological microscopy to excite fluorescence. However, most less-expensive biology microscopes such as those used in education and clinic applications are not equipped with epi-fluorescence illuminators. Even in more sophisticated biology microscopes EPI illuminators are only an expensive option. To make epi-fluorescence illumination available on more microscopes, an inexpensive illuminator that can fit on virtually any microscopes without alteration is called for.


Through-eyetube(TET) EPI illuminators are a unique class of microscope illumination solutions developed and patented by Mightex Systems. A TET-EPI illuminator sends light through one of the eyetubes of a typical binocular head. Light passes through microscope objective and illuminates the sample under observation, same as in a conventional EPI illumination setup. The TET-EPI illuminator fits into an eye tube just like a regular eyepiece. Image can be observed through the other eyepiece or the phototube. Alternatively, TE-EPI illuminator can be inserted into the phototube. In this configuration, images maybe observed through both eyepieces if the microscope is equipped with a beamsplitter between the phototube and the eyetubes.



PERFORMANCE SPECIFICATIONS

Parameters	Specification	Unit
Optical power density at exit pupil	>2	mW/cm ²
Excitation filter center wavelength	470	nm
Excitation filter bandwidth	40	nm
Eyepiece magnification	10x	-
Emission filter wavelength	510nm long-pass	-

SCHEMATIC DIAGRAM of TET-EPI



For more details about this product, please call 1-416-840-4991 or email sales@mightex.com.

FluoFlux Fluorescence LED Illuminators for Stereo-Microscopes

Stereo fluorescence has become a mainstream tool for developmental and cellular biology examinations. With fluorescence labeling, researchers can use stereo-microscopes to observe living organisms, screen samples, as well as to sort or dissect samples. Despite the growing needs, most stereomicroscopes are not equipped with fluorescence capabilities. Only high-end research stereomicroscopes have fluorescence excitation and imaging options. These options are also fairly expensive, often running into tens of thousands of dollars.

Mightex has developed and patented FluoFlux™, a LED-based stereo fluorescence solution that can be retrofit on most stereo microscopes. The solution includes a reflective excitation source and an emission filter fixture (note: filter not included) that can be mounted on lens barrel of the microscope.



PERFORMANCE SPECIFICATIONS

Table 1. LED Emitter Specifications

Code	Description	Peak Wavelength (nm)	Iop (mA)	Vop (V)	Radiant Flux (mW)
0309	3W UV	395 ~ 410	700	3.5	350
0301	3W Royal Blue	455	1000	3.9	450
0302	3W Blue	470	1000	3.9	480
0303	3W Cyan	505	1000	3.9	280
0304	3W Green	530	1000	3.9	135
0305	3W Amber	590	1400	3.0	210
0306	3W Red-Orange	617	1400	3.0	660
0307	3W Red	625	1400	3.0	630
0208	2.4W Near Infrared	850	1000	1.8~2.4	375

Table 2. Emission Holder

Holder Code	Inner diameter (mm)	Emission Filter diameter (mm)
66	66	36
60	60	36
55	55	36
50	50	36
45	45	36
00	No holder	NA

FilterReader Spectrophotometer

Ever wondered what's in your filter cube? FilterReader spectrophotometer will give you the answer in a few seconds.

Identifying and confirming transmission characteristics of filters and beam splitters used in fluorescence microscopy is often critical especially when troubleshooting a system or an experiment. To address this need Mightex has developed the FilterReader, a unique mini spectrophotometer specifically designed for measuring spectral transmission of excitation filters, dichroic beam splitters, and emission filters. A high-resolution CCD spectrograph coupled with a solid-state light source provides high dynamic range measurements in a compact package. A sample cart holds excitation filters and emission filters at vertical position, and beam splitters at a 45-degree angle.

FilterReader is a fully integrated spectrophotometer with a footprint of a book. All one needs is a PC to run the control software, which communicates with the FilterReader via a USB interface. FilterReader software guides users through an intuitive work flow which leads to the measurement results after a few clicks. Results can be displayed in spectral transmission or optical density. Measurements of filters and beam splitters can be overlaid on a single chart. Various data saving and chart printing functions are supported.

Unlike traditional spectrophotometers, FilterReader spectrophotometers are small enough to be easily stored when not in use. The solid-state light source in the FilterReader also ensures a trouble-free use for years without the need to change the light source.



PERFORMANCE SPECIFICATIONS

Parameters	Specifications		Unit
Model	FRD-350-AU	FRD-385-AU	
Wavelength Range	350 ~ 720	385 ~ 720	nm
Dynamic Range	>1000:1 (OD3)	>1000:1 (OD3)	
Filter Size	f25.4 or smaller	f25.4 or smaller	mm
Beam splitter Size	up to 27 x 36	up to 27 x 36	mm
PC Interface	USB2.0	USB2.0	
DC Power Supply	12V/1A (included)	12V/1A (included)	
Dimensions	41.5 x 108 x 48.5	41.5 x 108 x 48.5	mm
Weight	1,275	1,275	g

High-Resolution High-Stability CCD Spectrometers

Compact CCD spectrometers are widely used in process control, environment monitoring, and scientific research applications. Mightex HRS series compact CCD spectrometers features a high-resolution 100mm Czerny-Turner optical platform coupled with a Toshiba 3648-element CCD array. The optimized optical path yields both high spectral resolution and high light collection efficiency.



Wavelength and amplitude stability is often a critical requirement for many spectrometer applications. All optical components in the HRS series spectrometers are mounted directly on a single-piece base without using screws. A box enclosure structure further increases stiffness of the base. The proprietary mounting method ensures high stability over time and temperature.

A 16-bit DAC is used to convert the analog signal from the CCD array into a digital stream. The electronics hardware also includes trigger input and four programmable digital I/Os for interfacing with other equipment such as a light source. The spectrometer is controlled through a USB2.0 interface which also supplies all the electric power needed to operate the spectrometer.

Standard software package includes a full-featured PC software as well as a software development kit (SDK) for further software development.

PERFORMANCE SPECIFICATIONS

Parameters	Specifications				Unit
Model	HRS-UV1-025	HRS-VIS-025	HRS-NIR-025	HRS-BD1-025	-
Optical Platform	f/4, Czerny-Turner				-
Focal Length	100				mm
Wavelength Range	200~400	390~780	600~1000	300~1050	nm
Resolution*	0.25	0.4	0.5	0.9	nm
Order Sorting Filter	Long pass	Long pass	Long pass	Spatially Variable Filter	-
Entrance Slit	5, 10, 25, 50, 100, 200 or no slit				μm
Input Fiber Connector	SMA 905				-
Input Fiber NA	0.22				-
Detector	Toshiba TCD1304AP Linear CCD Array				-
Pixel Number	3648				-
Pixel Size	8 x 200				μm
Pixel Well Depth	100,000				electron
Signal-to-noise Ratio	1000:1 (at full scale)				-
A/D Resolution	16				bit
Integration Time	0.1 to 6,500				ms
Frame Rate	up to 138				fps
GPIO	4 programmable digital I/Os				-
Trigger Input	Yes				-
PC Interface	USB 2.0				-
Trigger/GPIO Interface	DIN8				-
PC Operating System	Windows 2000/XP/Vista/7				-
Power Consumption	300 (at 5V)				mA
Dimensions	138 x 108 x 37				mm
Weight	510				g

WARRANTY TERMS AND CONDITIONS

MIGHTEX's end user products and OEM modules are warranted to be free from defects in materials and workmanship for a period of 12 months and 6 months ("Warranty Period"), respectively, from the date of initial shipment. MIGHTEX'S liability under this warranty is limited to replacing any defective parts at its expense. MIGHTEX shall warrant the replacement products for the remainder of the original warranty period.

This warranty will not apply to those products: (i) repaired or altered other than in accordance with MIGHTEX's product specifications or written approval by MIGHTEX's duly authorized officer, or (ii) abused, misused, improperly handled in use or storage, or used in an unauthorized or improper manner or without following written procedures supplied by MIGHTEX, or (iii) original identification markings or labels have been removed, defaced or altered, or (iv) any other claims not arising directly from material defects in material or workmanship.

For board-level products (e.g. board-level cameras etc.) on which electronics components are exposed, one should be properly grounded when handling the products, in order to avoid Electro-Static Discharge (ESD) related damages. ESD related damages (which are considered as mis-handling) are not covered by product warranty.

Special contracts or contracts for nonstandard products may have modified terms of warranty and, in such cases, the terms as stated in the individual contract must be signed by the duly authorized officer of MIGHTEX and will supersede the standard terms.

MIGHTEX will make final determination as to cause or existence of defect and, at its option repair or replace the products, which prove to be defective during the warranty period. Products replaced under warranty will be warranted only for the balance of the warranty period from the original supplied equipment.

This warranty extends only to the original purchaser of the equipment from MIGHTEX. The purchaser must notify MIGHTEX within 15 days of first noticing the defect and promptly return the defective product upon receipt of RMA number(s) before expiration of the warranty period.

Products believed by purchaser to be defective shall be returned to MIGHTEX, with transportation and insurance prepaid by purchaser. Repaired or replaced products will be returned to purchaser by MIGHTEX, FOB city destination within the Continental United States and Canada, Transportation beyond these limits will be charged to purchaser.

The warranty set out in above paragraph is the exclusive warranty made by MIGHTEX and is in lieu of all other warranties (except for specific product performance warranties), whether written, oral, or implied, including any warranty of merchantability or fitness for a particular purpose, and shall be CUSTOMER'S sole remedy and MIGHTEX'S sole liability on contract or warranty of otherwise for the products. This warranty shall not be modified or amended without the written approval of an officer of MIGHTEX.

IN NO EVENT SHALL MIGHTEX BE LIABLE FOR ANY INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM THE FAILURE TO PERFORM UNDER THIS AGREEMENT OR THE FURNISHING, PERFORMANCE OR USE OF ANY GOODS OR SERVICES SOLD PURSUANT HERETO, WHETHER DUE TO BREACH OF CONTRACT, BREACH OF WARRANTY, NEGLIGENCE OR OTHERWISE.

CANCELLATION AND RETURN POLICY

After acceptance by the Seller, the Buyer's order is not subject to cancellation, change, reduction in amounts, or suspension or delay of shipment, except with the Seller's written consent. Cancelled orders consented to by the Seller are subject to cancellation charges, based upon all expenses incurred by the Seller up to the time of cancellation, including a reasonable profit on such expenses. This includes excess inventory of stock items that have been specifically allocated to the Buyer and exceed normal Seller inventory levels for a period. Seller may cancel the order in the following cases: the seller's costs have exceeded selling price to Buyer; the Seller is unable to provide product due to factors beyond the Seller's control.

If shipments are delayed by Buyer, payment shall become due on the date when Seller is prepared to make shipment unless otherwise agreed to in writing by Seller at time of incident. Products held for the Buyer shall be at the risk and expense of the Buyer.

To return Catalog Products, Buyer must contact Seller within 14 days of receipt to attain a Return Merchandise Authorization (RMA) Number. A 15% or a minimum of US\$50 restocking fee will be charged, and the Buyer is responsible for the shipping cost. After 14 days, the RMA number is still required; however, unless otherwise provided in a specific warranty, only products found to be defective are accepted for repair or replacement within the Warranty Period (as stated in Mightex's "Warranty Terms and Conditions") from invoice date. Return authorizations are valid 14 days from issue. All returned Products are subject to inspection and approval by Seller before issuing credit or replacement. Products displaying use, misuse, or other damage will not be accepted. Products must be in original manufacturer's packaging. Products should be packaged in a separate carton with the return authorization number clearly marked on this shipping carton. Refunds are issued by check for pre-paid merchandise. Credit order refunds are issued to the credit card or Seller's open account to which the Product was charged.

IMPORTANT:

Other than as expressly set forth above or as contained in any express warranties provided with products and the extent permitted by law, the seller makes no warranties, express or implied, including warranties of merchantability or fitness of the merchandise for any particular purpose. The seller shall not be liable for loss or damage resulting from the use or performance of the products. In no event shall the seller be liable to the buyer or its customers for any special, indirect, incidental, exemplary or punitive damages resulting from loss of use, interruption of business, or lost profits, even if the seller has been advised of the possibility of such damages. Notwithstanding the foregoing, in no event shall the seller's aggregate liability to the buyer and its customers exceed the amount paid for products purchased under this agreement in the prior 6 months.

GOVERNING LAW:

Sales transactions are governed by the internal law of the Province of Ontario, Canada. and you agree to the exclusive jurisdiction of the Ontario courts in the event of any dispute.

TAXES:

We collect sales tax on shipments to: (a) all Canadian provinces and territories, and (b) CA addresses, unless the purchaser supplies us with a signed official state resale certificate or sales tax exemption certificate with the order. Purchasers outside the aforementioned states/provinces are responsible for their state or local sales/use tax, if any, and any other taxes payable by reason of this transaction.

WARRANTY AND SAFETY INFORMATION:

Please refer to Mightex's Warranty Terms and Conditions.

REPAIRS:

If your unit is out of warranty but needs repairs, please write to us for repair information.

Mightex Systems

2343 Brimley Road, # 868

Toronto, Ontario M1S 3L6, Canada

Tel: 1-925-218-1885 or 1-416-840-4991

Fax: 1-416-840-6541

Email: sales@mightex.com

www.mightex.com or www.mightexsystems.com

This copy belongs to:

