


Doubly Curved Crystals

Highly monochromatic x-ray beams focused in three dimensions can be obtained by using doubly curved single crystals (DCCs). Crystals of mica, graphite, Si, Ge, and others are used for an x-ray energy range of ~1.5 to 22keV. Crystal optics reflect x-rays based on Bragg diffraction. The DCC accurately images micron-sized as well as large x-ray point sources.



TYPICAL FLUX AND SPOT SIZES

High-power rotating-anode sources:

- **Flux:** 1×10^{11} photons per second at the spot
- **Spot size:** ~300 μm x 200 μm

Microfocus sources:


- **Flux:** 1×10^9 photons per second at the spot
- **Spot size:** from 20 μm up to 150 μm (equivalent to source size)

FEATURES:

- 3-Dimensional point-to-point focusing of x rays
- Highly monochromatic beam
- Large capture angle
- High flux-density gain
- High angular uniformity of beam

BENEFITS:

- Extremely low background
- High detection sensitivity
- Detection limits in the ppb range
- Large working distance
- Rapid data collection
- Applications in MMXRF, XRR, and TXRF




MONOCHROMATIC MICRO X-RAY FLUORESCENCE ANALYSIS (MMXRF)

FIGURE 1

Comparison of a Mo DCC and a pinhole for elemental analysis of concentrated air particulates. Mo excitation (40kV, 20W, 600s)

The detection sensitivity provided by the DCC is greatly enhanced relative to a pinhole. The intensity of the pinhole spectrum was scaled from a 2 mm diameter aperture to a diameter equivalent to 50 μm. The pinhole and the DCC focal spot was located the same distance from the source (120 mm). A PIN detector with a 25 mm² active area was positioned 15 mm from the sample to collect both spectra.



— Pinhole, 50 micron diameter
 — DCC optic, 50 micron focal spot

MMXRF measurements below were done with an intense, monochromatic, focused beam, a Si(Li) detector, and a microfocus x-ray source (30kV, 3W). Results for high purity Ca and Ti specimens are shown in Table 1; results for mica are shown in Figure 2.

TABLE 1

Signal-to-background ratios for Ca and Ti samples Cu K α excitation (30kV, 3W, 500s).

Element	Signal*	Background*	S/B
Ca	1.8×10^7	140	1.3×10^5
Ti	3×10^7	200	1.5×10^5

The signal-to-background ratios for Ca and Ti measurements are more than two orders of magnitude higher than those provided by standard XRF techniques.

* Signal and background counts are integrated in 10 channels (20eV/channel)

FIGURE 2

Fluorescence spectrum of a mica sample using Cu K α excitation (30kV, 3W, 500s).

Concentrations at the low ppm level were detected in 500s with the source operated at only 3W. Under these conditions the primary beam still delivered a total of $\sim 6 \times 10^7$ photons.





FIGURE 3**Comparison of spectra from a polycapillary optic and Cu DCC, both with 50µm focal spot size.**

Si DCC–30kV, 20W, 400s

Polycapillary optic–20kV, 20W, 120s.

Spectra collected from a Lucite® plastic disk. The data shows the background suppression benefit of using DCC optics.



Low background and low detection-limits can be achieved with doubly curved crystals.

TABLE 2**Examples of different DCCs integrated with microfocus x-ray sources:**


Optic	Energy (keV)	Source Power (W)	θ Bragg	Capture Angle (sr)	Nominal focal spot size, FWHM (μm)	Flux (cps)	Working Distance (mm)
Cr1	5.4	14	35.1	0.03	80	2×10^9	120
Cu1	8.0	14	14.2	0.015	50	1×10^9	150
Cu2	8.0	50	14.2	0.01	150	2×10^9	150
Cu3	8.0	14	22.8	0.01	50	3×10^8	135
W1	8.4	10	22.6	0.01	20	1×10^8	80
Mo1	17.5	14	10.6	0.01	60	1×10^8	120

MONOCHROMATIC CONVERGENT X-BEAM®

For x-ray analysis requiring high flux-density gain and micron-sized x-ray focal spots, XOS® offers DCC optics integrated with microfocus x-ray sources. For example, we offer the following nominal specification:

Flux (Mo $K\alpha$): 3.0×10^8 cps at 50 WSpot (FWHM): 200 μm

Please contact the sales department for information.





Thin Film/Surface Analysis

X-RAY REFLECTOMETRY (XRR)

XRR data for 800Å TiN film on Si wafer using Tungsten L α line

- 3.5W (35kV, 0.1mA) source setting; 100s data collection
- The data was collected using a PSD by scanning a 65 μ m slit positioned 400 mm from the sample. A remarkable angular resolution of $\sim 0.01^\circ$ and very high count rate was generated using a low-power Tungsten x-ray source with a 40 μ m spot size.





Schematic diagram of XRR geometry.

TOTAL REFLECTION X-RAY FLUORESCENCE (TXRF)

Spectrum of a residue on the surface of a silicon wafer, showing 3pg Br detection using Mo K α excitation

- 40W (40kV, 1mA) source setting; 500s data collection
- The plot shows a low femto-gram detection limit of about 2×10^8 atoms with a compact, low-power x-ray source. A small beam size of 100 μ m is achieved using a Si Doubly Curved Crystal for particle analysis and mapping.



Schematic diagram of TXRF geometry.

X-RAY
OPTICAL SYSTEMS INC.

The World Leader in X-Ray Optics and Beams

XOS® offers standard and custom designs for your analytical needs.



光技術をサポートする
株式会社オプトサイエンス

<http://www.optoscience.com>

東京本社 〒160-0014 東京都新宿区内藤町1番地 内藤町ビルディング TEL:03-3356-1064
大阪営業所 〒532-0011 大阪市淀川区西中島7-7-2 新大阪ビル西館 TEL:06-6305-2064
名古屋営業所 〒450-0002 名古屋市中村区名駅2-37-21 東海ソフトビル TEL:052-569-6064

E-mail : info@optoscience.com