# Andor iXon Ultra Blue

OXFORD  
INSTRUMENTS

ANDOR

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## Advanced Performance Specifications ※2

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※2 Figures are typical unless otherwise stated.

※4 The dark current measurement is averaged over the sensor area excluding any regions of blemishes.

※5 Using Electron Multiplication the iXon is capable of detecting single photons, therefore the true camera detection limit is set by the number of 'dark' background events. These events consist of both residual thermally generated electrons and Clock Induced Charge (CIC) electrons (also referred to as Spurious Noise), each appearing as random single spikes above the read noise floor. A thresholding scheme is employed to count these single electron events and is quoted as a probability of an event per pixel. Acquisition conditions are full resolution and max frame rate (30 MHz readout; frame transfer mode; 1.1  $\mu\text{s}$  vertical clock speed; x 1000 EM gain; 10 ms exposure; -95°C).

※6 The EM register on CCD201 sensors has a linear response up to ~400,000 electrons and a full well depth of ~730,000 electrons.

※7 Readout noise is for the entire system. It is a combination of sensor readout noise and A/D noise. Measurement is for Single Pixel readout with the sensor at a temperature of -75°C and minimum exposure time under dark conditions. Under Electron Multiplying conditions, the effective system readout noise is reduced to sub 1 e- levels.

※8 Linearity is measured from a plot of counts vs. exposure time under constant photon flux up to the saturation point of the system, at 10 MHz readout speed.

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標準モード ※3,9								
		1018 x 1018	506 x 506	256 x 256	128 x 128	1018 x 100	1018 x 32	1018 x 1
Binning	1 x 1	26	50	95	171	220	498	1163
	2 x 2	50	92	170	185	368	699	-
	4 x 4	92	167	281	426	552	870	-

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
※3 At 30 MHz overclocked pixel readout rate, thermal dissipation from the sensor is higher since a greater proportion of time is spent vertical shifting, and it is necessary to set a higher sensor cooling temperature at this rate. Furthermore, stable cooling performance will depend on other variables such as vertical clock speed, Region of Interest size (Standard or Crop Mode) and ambient temp. As such, user testing is advised to determine the stable sensor cooling temperature for selected conditions. Status of temperature stability is apparent through the acquisition software.

※9 All measurements are made at 30 MHz pixel readout speed with 0.6 μs vertical clock speed. It also assumes internal trigger mode of operation. Standard and Crop Mode frame rates shown are for 'Corner Tethered' ROIs, with 'Optically Centred' ROI frame rates shown within brackets.

※10 All measurements are made at 17 MHz pixel readout speed with 0.3 μs vertical clock speed. It also assumes internal trigger mode of operation. Standard and Crop Mode frame rates shown are for 'Corner Tethered' ROIs, with 'Optically Centred' ROI frame rates shown within brackets.

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
Creating The Optimum Product for You

## Step 1. Choose the camera type

Camera Type

Description	Code
iXon Ultra 888: 1018 x 1018 EMCCD, max. 30 MHz, with USB 3.0	<b>DU-888U3-CS0-NBB</b>
iXon Ultra 897: 506 x 506 EMCCD, max. 17 MHz, with USB 2.0	<b>DU-897U-CS0-NBB</b>


## Step 2. Select an alternative camera window (optional)

Camera Window

The standard window has been selected to satisfy most applications. However, other options are available. The alternative camera window code must be specified at time of ordering.  
To view and select other window options please refer to the '[Camera Windows Supplementary Specification Sheet](#)' which gives the transmission characteristics, product codes and procedure for entering the order. Further detailed information on the windows is in the Technical note – '[Camera Windows: Optimizing for Different Spectral Regions](#)'.

## Step 3. Select the required accessories

Accessories

Description	Order Code	Description	Order Code
SRRF-Stream Dell Workstation (English), pre-installed with a recommended and tested GPU card, alongside SRRF-Stream enabled MicroManager and Andor SDK2 with SRRF-Stream.	<b>WKST-SRRF-9ZY</b>	Re-circulator for enhanced cooling performance	<b>XW-RECR</b>
Monitor (optional) - Dell UltraSharp U3417W - 34.14" Curved LED	<b>FUS-MNTR-34W</b>	Oasis 160 Ultra compact chiller unit (tubing to be ordered separately)	<b>ACC-XW-CHIL-160</b>
Dell UltraSharp UP3017 - 30" with PremierColor	<b>FUS-MNTR-30</b>	6 mm tubing options for ACC-XW-CHIL-160 (2x2.5 m or 2x5m lengths)	<b>ACC-6MM-TUBING-2X2.5/ ACC-6MM-TUBING-2X5M</b>
OptoMask accessory, used to mask unwanted sensor area during Crop Mode acquisition (refer to OptoMask Specification Sheet for further information).	<b>OPTMSK-L/ OPTMSK-OC-L/ OPTMSK-OC-S</b>	C-mount to Nikon F-mount adapter	<b>OA-CNAF</b>
		C-mount to Olympus adapter	<b>OA-COFM</b>
		C-mount to T-mount adapter	<b>OA-CTOT</b>
		15 m Active USB 3.0 connector cable (power supply not required) Icron for Ultra 888	<b>ACC-ASE-06887</b>
		50 m Fibre Optic USB 3.0 extender solution inc. power supply (Adnaco) for Ultra 888	<b>ACC-ASE-08762</b>
		100 m Fibre Optic USB 3.0 extender solution inc. power supply (Adnaco) for Ultra 888	<b>ACC-ASE-07860</b>

## Step 4. Select the required software


Software

The iXon Ultra series requires one of the following software options:  
**Solis Imaging:** A 32-bit and fully 64-bit enabled application for Windows (8, 8.1, 10 and 11) offering rich functionality for data acquisition and processing. AndorBasic provides macro language control of data acquisition, processing, display and export.  
**Andor SDK:** A software development kit that allows you to control the Andor range of cameras from your own application. Available as 32 and 64-bit libraries for Windows (8, 8.1, 10 and 11), compatible with C/C++, C#, Delphi, VB.NET, LabVIEW and Matlab. Linux SDK compatible with C/C++.  
**Third party software compatibility,** drivers are available for a variety of [third party imaging packages](#).

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


## iXon Ultra 888

## iXon Ultra 897

※1 Assembled in a state-of-the-art cleanroom facility, Andor's UltraVac™ vacuum process combines a permanent hermetic vacuum seal (no o-rings), with a stringent protocol to minimize outgassing, including use of proprietary materials.

※12 iXon Ultra 888 should work with any modern USB 3.0 enabled PC/laptop, as every USB 3.0 port should have its own host controller. iXon Ultra 888 also ships with a USB 3.0 PCI card as a means to add a USB 3.0 port to an older PC, or as a diagnostic aid to interoperability issues.

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#### Items shipped with your iXon Ultra 888:

1x Andor ACZ-03463: 2m Multi I/O timing cable, offering Fire, External Trigger, Shutter and Arm  
1x 3m USB 3.0 cable Type A to Type B  
PCIe USB 3.0 Card Adapter (2-Port) ※12  
1x Power supply unit with mains cable  
1x Quick Start guide  
1x Electronic copy of user manuals  
1x SRRF-Stream Quick Start guide (if applicable)  
1x Individual system performance booklet

#### Items shipped with your iXon Ultra 897:

1x Andor ACZ-03463: 2 m Multi I/O timing cable, offering Fire, External Trigger, Shutter and Arm  
1x 3m USB 2.0 cable Type A to Type B  
1x Power supply unit with mains cable  
1x Quick Start guide  
1x Electronic copy of user manuals  
1x SRRF-Stream Quick Start guide (if applicable)  
1x Individual system performance booklet

#### Recommended Computer Requirements:

- 3.0 GHz single core or 2.6 GHz multi core processor
- 2 GB RAM
- 100 MB free disc space to install software (at least 1 GB recommended for data spooling)
- USB 3.0 Super Speed Host Controller capable of a sustained rate of 60MB/s for iXon Ultra 888
- USB 2.0 High Speed Host Controller capable of sustained rate of 40MB/s for iXon Ultra 897
- Solid-state drive (SSD) capable of a minimum sustained write speed of 100MB/S for spooling data
- Windows (8.1, 10 and 11) or Linux
- SRRF-Stream+ - If selected, the PC requires a Nvidia GPU card. See page 10 for further details.

#### Operating & Storage Conditions

- Operating Temperature: -20°C to 30°C ambient
- Relative Humidity: < 70% (noncondensing)
- Storage Temperature: -25°C to 50°C

#### Power Requirements

- Please refer to page 10