




FEMTO® PRODUCT OVERVIEW

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

FEMTO®



光技術をサポートする
株式会社オプトサイエンス
<http://www.optoscience.com>

東京本社 〒160-0014 東京都新宿区内藤町1番地 内藤町ビルディング
 TEL: 03 (3356) 1064 FAX: 03 (3356) 3466 E-mail: info@optoscience.com
 大阪支店 〒532-0011 大阪市淀川区西中島7-7-2 新大阪ビル西館
 TEL: 06 (6305) 2064 FAX: 06 (6305) 1030 E-mail: osk@optoscience.com
 名古屋営業所 〒450-0002 名古屋市中村区名駅2-37-21 東海ソフトビル
 TEL: 052 (569) 6064 FAX: 052 (569) 8064 E-mail: ngo@optoscience.com

FEMTO AMPLIFIER SELECTION GUIDE




TABLE OF CONTENTS

| | |
|---|----|
| High Speed Charge Amplifier ■ Model HQA-15M-10T | 4 |
| Ultra Low Noise Current Amplifiers ■ Series LCA | 5 |
| High Speed Current Amplifiers ■ Series HCA | 5 |
| Sub Femto Ampere Current Amplifier ■ Model DDPCA-300 | 6 |
| Variable Gain Low Noise Current Amplifier ■ Model DLPCA-200 | 6 |
| Variable Gain High Speed Current Amplifier ■ Model DHPCA-100 | 7 |
| Variable Gain 100 kHz Voltage Amplifiers ■ Series DLPVA-100 | 8 |
| Wideband Voltage Amplifiers ■ Series HVA | 8 |
| Variable Gain 100/200 MHz Voltage Amplifiers ■ Series DHPVA | 9 |
| Logarithmic Wideband Voltage Amplifier ■ Model HLVA-100 | 9 |
| High Speed GHz Amplifiers ■ Series HSA-X | 10 |
| High Speed GHz Amplifiers ■ Series HSA-Y | 10 |
| Variable Gain GHz Voltage Amplifiers ■ Series DUPVA | 11 |
| Miniature Lock-In Amplifier Modules ■ Series LIA-MV-150 | 12 |
| Lock-In Amplifier Modules ■ Series LIA-MV(D)-200 | 12 |
| Single-Board Lock-In Amplifiers ■ Series LIA-BV(D)-150 | 13 |
| Multi-Channel Lock-In Amplifier Rack ■ Series SC-LIA-S | 13 |
| Femtowatt Photoreceiver ■ Series FWPR-20 | 15 |
| 400 kHz Low Noise Photoreceiver ■ Series LCA-S-400K | 15 |
| 500 kHz Variable Gain Photoreceiver ■ Series OE-200 | 16 |
| 200 MHz High Speed Photoreceiver ■ Series HCA-S-200M | 16 |
| 400 MHz High Speed Photoreceiver ■ Series HCA-S-400M | 17 |
| 2 GHz High Speed Photoreceiver ■ Series HSA-X-S | 17 |
| Power Supply ■ Series PS-15 | 18 |
| USB Control Interface ■ Model LUCI-10 | 18 |
| Custom-Designed Solutions | 19 |

High Speed Charge Amplifier ■ Model HQA-15M-10T



Model HQA-15M-10T

- High Gain of 10 V/pC
- Wide Operating Range from 250 Hz to 15 MHz
- Low Input Noise of 40×10^{-21} C/√Hz and 700 pV/√Hz
- Ideal for AC Coupled Charge Sources Like Pyro- and Piezoelectric Detectors, Tuning Fork Quartz Crystals and Length Extension Resonators
- Typical Applications are Atomic Force Microscopy, Optical Measurements and Charged Particle Beam Monitoring

| Model | HQA-15M-10T |
|--------------------------------|---|
| Charge Gain | 10^{13} V/C |
| Equivalent Current Gain | 1.6×10^6 V/A @ 1 MHz Sinusoidal Input Signal |
| Lower Cut-Off Frequency | 250 Hz (AC only) |
| Upper Cut-Off Frequency | 15 MHz |
| Input Charge Noise | 40×10^{-21} C/√Hz |
| Equivalent Input Current Noise | 250 fA/√Hz @ 1 MHz Sinusoidal Input Signal |
| Input Voltage Noise | 700 pV/√Hz @ 1 MHz |
| Input Impedance | $1 \text{ G}\Omega // 1 \text{ nF}$ |
| Effective AC Input Impedance | $20 \Omega @ 1 \text{ MHz Sinusoidal Input Signal}$ |
| Output Performance | 10 Vp-p @ $\geq 1 \text{ M}\Omega$ Load |
| Power Supply | $\pm 15 \text{ V}, \pm 35 \text{ mA Typ.}$ |
| Case | 115 x 60 x 45 mm (L x W x H), Weight 200 g (0.44 lbs) |

Maximum recommended source capacitance is 1 nF. Not suited for sources producing an average DC background current. Output short-circuit protected. Power supply via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.

Ultra Low Noise Current Amplifiers ■ Series LCA




Series LCA

- Ultra Low Input Noise down to 0.2 fA/√Hz
- Gain up to 10^{13} V/A
- Bandwidth DC up to 400 kHz
- Bandwidth and Gain Independent of Source Capacitance up to 10 nF (1 nF for LCA-400K-10M)
- Compact and Highly Shielded Case for Use Close to the Signal Source
- Ideal for Photodiodes, STMs and Ionization Detectors

| Model | - 3 dB Bandwidth (DC ...) | Noise Current | Transimpedance (Gain) | Rise/Fall Time |
|--------------|---------------------------|---------------|---------------------------|----------------|
| LCA-2-10T | 2 Hz | 0.2 fA/√Hz | 10^{12} & 10^{13} V/A | 200 ms |
| LCA-30-1T | 30 Hz | 0.5 fA/√Hz | 1×10^{12} V/A | 12 ms |
| LCA-30-200G | 30 Hz | 0.5 fA/√Hz | 2×10^{11} V/A | 12 ms |
| LCA-200-100G | 200 Hz | 1.5 fA/√Hz | 1×10^{11} V/A | 2 ms |
| LCA-200-10G | 200 Hz | 1.5 fA/√Hz | 1×10^{10} V/A | 2 ms |
| LCA-1K-5G | 1 kHz | 3 fA/√Hz | 5×10^9 V/A | 400 μs |
| LCA-2K-2G | 2 kHz | 4.5 fA/√Hz | 2×10^9 V/A | 200 μs |
| LCA-4K-1G | 4 kHz | 6.5 fA/√Hz | 1×10^9 V/A | 100 μs |
| LCA-10K-500M | 10 kHz | 10 fA/√Hz | 5×10^8 V/A | 40 μs |
| LCA-20K-200M | 20 kHz | 14 fA/√Hz | 2×10^8 V/A | 20 μs |
| LCA-40K-100M | 40 kHz | 19 fA/√Hz | 1×10^8 V/A | 10 μs |
| LCA-100K-50M | 100 kHz | 30 fA/√Hz | 5×10^7 V/A | 4 μs |
| LCA-200K-20M | 200 kHz | 40 fA/√Hz | 2×10^7 V/A | 2 μs |
| LCA-400K-10M | 400 kHz | 65 fA/√Hz | 1×10^7 V/A | 1 μs |

Bandwidth and rise/fall times are independent of detector capacitance. Guaranteed and tested up to 10 nF for each amplifier (1 nF for LCA-400K-10M). Output voltage ± 10 V @ ≥ 1 M Ω load. Offset adjustable by trimpot. Output short-circuit protected. Power supply ± 15 V via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.

High Speed Current Amplifiers ■ Series HCA




Series HCA

- High Bandwidth from DC up to 400 MHz
- Optimized Models for High Source Capacitance up to 2 nF
- Adjustable Bias Voltage and Offset
- Ideal for Fast and Large Area Photodiodes, PMTs and Ionization Detectors

| Model | - 3 dB Bandwidth (DC ...) | Noise Current | Transimpedance (Gain) | Rise/Fall Time | Max. Detector Capacitance |
|----------------|---------------------------|---------------|-----------------------|----------------|---------------------------|
| HCA-1M-1M | 1 MHz | 270 fA/√Hz | 1×10^6 V/A | 350 ns | 50 pF |
| HCA-1M-1M-C | 1 MHz | 3.5 pA/√Hz | 1×10^6 V/A | 350 ns | 2 nF |
| HCA-2M-1M | 2 MHz | 340 fA/√Hz | 1×10^6 V/A | 180 ns | 25 pF |
| HCA-2M-1M-C | 2 MHz | 3.5 pA/√Hz | 1×10^6 V/A | 180 ns | 1 nF |
| HCA-4M-500K | 4 MHz | 490 fA/√Hz | 5×10^5 V/A | 90 ns | 15 pF |
| HCA-4M-500K-C | 4 MHz | 3.5 pA/√Hz | 5×10^5 V/A | 90 ns | 500 pF |
| HCA-10M-100K | 10 MHz | 1.1 pA/√Hz | 1×10^5 V/A | 35 ns | 15 pF |
| HCA-10M-100K-C | 10 MHz | 3.5 pA/√Hz | 1×10^5 V/A | 35 ns | 150 pF |
| HCA-20M-100K-C | 20 MHz | 3.5 pA/√Hz | 1×10^5 V/A | 18 ns | 50 pF |
| HCA-40M-100K-C | 40 MHz | 3.7 pA/√Hz | 1×10^5 V/A | 10 ns | 30 pF |
| HCA-100M-50K-C | 100 MHz | 3.8 pA/√Hz | 5×10^4 V/A | 3.5 ns | 20 pF |
| HCA-200M-20K-C | 200 MHz | 4.9 pA/√Hz | 2×10^4 V/A | 1.9 ns | 8 pF |
| HCA-400M-5K-C | 400 MHz | 21 pA/√Hz | 5×10^3 V/A | 1.0 ns | 10 pF |

The maximum detector capacitance listed above means that up to this value the specified bandwidth ($\pm 15\%$) is guaranteed. Larger capacitances are also possible but will influence the bandwidth. Output voltage ± 1.5 V @ 50 Ω load. Offset adjustable by trimpot. Output short-circuit protected. Adjustable bias output (-12 V ... +12 V) for biasing photodetectors. Power supply ± 15 V via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.

Variable Gain Sub Femto Ampere Current Amplifier ■ Model DDPKA-300



Model DDPKA-300


- 0.4 fA Peak-Peak Noise
- Variable Gain over 10 Decades from 10^4 to 10^{13} V/A for Sub-fA to mA Measurements
- Sub-fA Drift
- Optimized for DC to 400 Hz Measurements
- Adjustable Bias Voltage on Input for DUT Biasing
- Manual and Remote Control Operation
- Ideal for I/V Characterization of MOS and JFET Structures, Ultra Low Current and High Resistance Measurements, Quantum and Biotech Experiments, and as Easy to Use Sub Femto Amp Add-on to an Existing Digital Voltmeter or A/D Converter

| Model | DDPKA-300 | | | | | | | | | |
|--------------------------|--|--------|---------|---------|--------|--------|-----------|-----------|-----------|-----------|
| Transimpedance [V/A] | 10^4 | 10^5 | 10^6 | 10^7 | 10^8 | 10^9 | 10^{10} | 10^{11} | 10^{12} | 10^{13} |
| Bandwidth (-3 dB) [Hz] | 400 | 400 | 400 | 400 | 150 | 150 | 20 | 20 | 1 | 1 |
| Rise Time (10% - 90 %) | 0.8 ms | 0.8 ms | 0.8 ms | 0.8 ms | 2.3 ms | 2.3 ms | 17 ms | 17 ms | 350 ms | 350 ms |
| Equ. Input Noise [fV/Hz] | 45 pA | 45 pA | 0.45 pA | 0.45 pA | 15 fA | 15 fA | 1.3 fA | 1.3 fA | 0.2 fA | 0.2 fA |
| Low Pass Filter | Switchable to 0.1 Hz, 0.7 Hz or Full Bandwidth | | | | | | | | | |
| Output Performance | ± 10 V @ ≥ 1 M Ω Load | | | | | | | | | |
| Bias Voltage | ± 10 V, Adjustable by Trimpot or External Control Voltage, max. 10 mA, Connected to Center Pin of BNC Input Socket | | | | | | | | | |
| Power Supply | ± 15 V, + 70 mA / - 15 mA Typ. | | | | | | | | | |
| Control Interface | 4 Opto-Isolated Digital Inputs, TTL/CMOS Compatible, Analog Bias Control Voltage Input | | | | | | | | | |
| Case | 170 x 60 x 45 mm (L x W x H), Weight 320 g (0.74 lbs) | | | | | | | | | |

Offset adjustable by trimpot. LED overload indication. Input protected against ± 2 kV transients. Output short-circuit protected. Power supply via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.

Variable Gain Low Noise Current Amplifier ■ Model DLPCA-200

- Variable Gain over 9 Decades from 10^3 to 10^{11} V/A
- Low Input Noise down to 4 fA/Hz
- Bandwidth up to 500 kHz
- Adjustable Offset and Bias Voltage
- Switchable AC/DC Coupling and 10 Hz Low Pass Filter for Precise DC Measurements
- Manual and Remote Control Operation
- Designed for Photodetectors, STMs, Ionization Detectors, as Pre-Amplifier for Lock-Ins, A/D Converters and for General Lab Use




Model DLPCA-200

| Model | DLPCA-200 | | | | | | | | | | | | | |
|--------------------------|---|--------|--------|-------------|-----------|------------|-------------|------------|--------|--------|-------------|-----------|------------|-------------|
| Performance Range | Low Noise | | | | | | | High Speed | | | | | | |
| Transimpedance [V/A] | 10^3 | 10^4 | 10^5 | 10^6 | 10^7 | 10^8 | 10^9 | 10^5 | 10^6 | 10^7 | 10^8 | 10^9 | 10^{10} | 10^{11} |
| Bandwidth (-3 dB) [kHz] | 500 | 500 | 400 | 200 | 50 | 7 | 1.1 | 500 | 500 | 400 | 200 | 50 | 7 | 1.1 |
| Rise Time (10% - 90 %) | 700 ns | 700 ns | 900 ns | 1.8 μ s | 7 μ s | 50 μ s | 300 μ s | 700 ns | 700 ns | 900 ns | 1.8 μ s | 7 μ s | 50 μ s | 300 μ s |
| Equ. Input Noise [fV/Hz] | 20 pA | 2.3 pA | 450 fA | 130 fA | 43 fA | 13 fA | 4.3 fA | 13 pA | 1.8 pA | 440 fA | 130 fA | 43 fA | 13 fA | 4.3 fA |
| Low Pass Filter | Switchable to 10 Hz | | | | | | | | | | | | | |
| Output Performance | ± 10 V @ ≥ 1 M Ω Load | | | | | | | | | | | | | |
| Bias Voltage | Adjustable ± 10 V, max. 22 mA, Connected to Shield of BNC Input Socket, Switchable to GND | | | | | | | | | | | | | |
| Power Supply | ± 15 V, + 120 mA / - 80 mA Typ. | | | | | | | | | | | | | |
| Control Interface | 5 Opto-Isolated Digital Inputs, TTL/CMOS Compatible, Analog Offset Control Voltage Input | | | | | | | | | | | | | |
| Case | 170 x 60 x 45 mm (L x W x H), Weight 320 g (0.74 lbs) | | | | | | | | | | | | | |

Bandwidth and rise/fall times are independent of detector capacitance up to 1 nF. Offset adjustable by trimpot or external control voltage. LED overload indication. Input protected against ± 3 kV transients. Output short-circuit protected. Power supply via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.

Variable Gain High Speed Current Amplifier ■ Model DHPCA-100

- Variable Gain over 7 Decades from 10^2 to 10^8 V/A
- Bandwidth up to 200 MHz
- 1.8 ns Minimum Rise Time
- Switchable Signal Filters and AC/DC Coupling
- Adjustable Offset and Bias Voltage
- Manual and Remote Control Operation
- Ideal for Fast Photodetectors, as Pre-Amplifier for Oscilloscopes, RF Lock-Ins and for General Lab Use



Model DHPCA-100

| Model | DHPCA-100 | | | | | | | | | | | |
|---|---|--------|--------|-------------|-------------|-------------|------------|--------|--------|-------------|-------------|-------------|
| Performance Range | Low Noise | | | | | | High Speed | | | | | |
| Transimpedance [V/A] | 10^2 | 10^3 | 10^4 | 10^5 | 10^6 | 10^7 | 10^3 | 10^4 | 10^5 | 10^6 | 10^7 | 10^8 |
| Bandwidth (-3 dB) [MHz] | 200 | 80 | 14 | 3.5 | 1.8 | 0.22 | 175 | 80 | 14 | 3.5 | 1.8 | 0.22 |
| Rise Time (10% - 90 %) | 1.8 ns | 4.4 ns | 25 ns | 0.1 μ s | 0.2 μ s | 1.6 μ s | 2.0 ns | 4.4 ns | 25 ns | 0.1 μ s | 0.2 μ s | 1.6 μ s |
| Equ. Input Noise [$\sqrt{\text{V/Hz}}$] | 200 pA | 16 pA | 2.1 pA | 500 fA | 170 fA | 60 fA | 140 pA | 6.0 pA | 1.5 pA | 450 fA | 150 fA | 55 fA |
| Low Pass Filter | Switchable to 1 MHz, 10 MHz or Full Bandwidth | | | | | | | | | | | |
| Output Performance | ± 1 V @ 50 Ω Load | | | | | | | | | | | |
| Bias Voltage | Adjustable ± 10 V, max. 22 mA, Connected to Shield of BNC Input Socket, Switchable to GND | | | | | | | | | | | |
| Power Supply | ± 15 V, + 110 mA / - 90 mA Typ. | | | | | | | | | | | |
| Control Interface | 7 Opto-Isolated Digital Inputs, TTL/CMOS Compatible, Analog Offset Control Voltage Input | | | | | | | | | | | |
| Case | 170 x 60 x 45 mm (L x W x H), Weight 320 g (0.74 lbs) | | | | | | | | | | | |

Offset adjustable by trimpot or external control voltage. LED overload indication. Input protected against ± 3 kV transients. Output short-circuit protected. Power supply via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.

Variable Gain Current Amplifiers ■ Comparison

| Model | DDPCA-300 | DLPCA-200 | DHPCA-100 |
|---------------------------------|--|-------------------------------|-------------------------------|
| Gain Range | $10^4 - 10^{13}$ V/A | $10^3 - 10^{11}$ V/A | $10^2 - 10^8$ V/A |
| Max. Bandwidth | DC to 400 Hz | DC to 500 kHz | DC to 200 MHz |
| Min. Rise Time | 800 μ s | 700 ns | 1.8 ns |
| Min. Spectral Noise | 0.2 fA/ $\sqrt{\text{Hz}}$ | 4 fA/ $\sqrt{\text{Hz}}$ | 55 fA/ $\sqrt{\text{Hz}}$ |
| Min. Integrated Broadband Noise | 0.4 fA peak-peak or 0.06 fA rms | 200 fA peak-peak or 30 fA rms | 500 pA peak-peak or 75 pA rms |
| Main Differentiator | Sub Femto Ampere Sensitivity | Broad Application Range | MHz Speed |
| Common Features | Manual and Remote Control Operation Adjustable Offset and Bias Voltage Compact Design for Use Close to the Signal Source Extremely Well Shielded Housing to Avoid Noise Pick-Up | | |

Variable Gain 100 kHz Voltage Amplifiers ■ Series DLPVA-100

- Low Input Noise down to 400 pV/√Hz
- Variable Gain up to 100 dB (x 100,000)
- Single Ended or Differential Input with Common Mode Rejection (CMRR) up to 120 dB
- High Input Impedance up to 1 TΩ
- For Amplification of Low Frequency Signals from Low to Medium Impedance Sources




Model DLPVA-100-B-S

| Model | DLPVA-100-BUN-S | DLPVA-100-BLN-S | DLPVA-100-B-S | DLPVA-100-B-D | DLPVA 100-F-S | DLPVA-100-F-D |
|-------------------------|-------------------------|-------------------------|-------------------------|----------------------------|-------------------------|-------------------------|
| Input stage | Single Ended, Bipolar | Single Ended, Bipolar | Single Ended, Bipolar | True Differential, Bipolar | Single Ended, FET | True Differential, FET |
| Input | 1 kΩ, BNC | 1 MΩ, BNC | 1 MΩ, BNC | 1 MΩ, LEMO | 1 TΩ, BNC | 1 TΩ, LEMO |
| Lower Cut-Off Frequency | 1.5 Hz (AC only) | DC/1.5 Hz | DC/1.5 Hz | DC/1.5 Hz | DC/1.5 Hz | DC/1.5 Hz |
| Upper Cut-Off Frequency | 1/100 kHz | 1/100 kHz | 1/100 kHz | 1/100 kHz | 1/100 kHz | 1/100 kHz |
| Gain | 40/60/80/100 dB | 40/60/80/100 dB | 20/40/60/80 dB | 20/40/60/80 dB | 20/40/60/80 dB | 20/40/60/80 dB |
| Input Noise Voltage | 0.4 nV/√Hz | 0.7 nV/√Hz | 2.4 nV/√Hz | 3.6 nV/√Hz | 5.5 nV/√Hz | 6.9 nV/√Hz |
| Input Voltage Drift | – | 0.5 μV/°C | 0.7 μV/°C | 0.7 μV/°C | 1.3 μV/°C | 1.3 μV/°C |
| CMRR | – | – | – | 120 dB Max. | – | 120 dB Max. |
| Output | 50 Ω, BNC | 50 Ω, BNC | 50 Ω, BNC | 50 Ω, BNC | 50 Ω, BNC | 50 Ω, BNC |
| Output Voltage | ± 10 V @ ≥ 1 MΩ | ± 10 V @ ≥ 1 MΩ | ± 10 V @ ≥ 1 MΩ | ± 10 V @ ≥ 1 MΩ | ± 10 V @ ≥ 1 MΩ | ± 10 V @ ≥ 1 MΩ |
| Control Interface | TTL/CMOS, Opto-Isolated | TTL/CMOS, Opto-Isolated | TTL/CMOS, Opto-Isolated | TTL/CMOS, Opto-Isolated | TTL/CMOS, Opto-Isolated | TTL/CMOS, Opto-Isolated |

Offset adjustable by trimpot or external analog control voltage. LED gain setting and overload indication. Output short-circuit protected. Power supply ± 15 V via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.

Wideband Voltage Amplifiers ■ Series HVA



Model HVA-200M-40-B

- Bandwidth DC up to 500 MHz
- Switchable Gain up to 60 dB (x 1,000)
- Very Low Noise down to 0.9 nV/√Hz
- True DC Coupling Allows Accurate Amplification of Digital Codes and Transients with Long Decay Time (No Baseline Shift or Signal Over/Under-Shoots)
- Ideal as Pre-Amplifier for Oscilloscopes and Transient Recorders, and for Time-Resolved Pulse and Transient Measurements

| Model | HVA-10M-60-B | HVA-10M-60-F | HVA-200M-40-B | HVA-200M-40-F | HVA-500M-20-B |
|-------------------------|----------------|----------------|---------------|---------------|---------------|
| Lower Cut-Off Frequency | DC/1 kHz | DC/1 Hz | DC/1 kHz | DC/1 Hz | DC |
| Upper Cut-Off Frequency | 10 MHz | 10 MHz | 200 MHz | 200 MHz | 500 MHz |
| Input Impedance | 50 Ω, BNC | 1 MΩ, BNC | 50 Ω, BNC | 1 MΩ, BNC | 50 Ω, BNC |
| Gain | 40/60 dB | 40/60 dB | 20/40 dB | 20/40 dB | 20 dB |
| Input Noise Voltage | 0.9 nV/√Hz | 4.7 nV/√Hz | 1.2 nV/√Hz | 4.5 nV/√Hz | 3.0 nV/√Hz |
| Input Voltage Drift | 1 μV/°C | 2 μV/°C | 1 μV/°C | 5 μV/°C | 10 μV/°C |
| Output Voltage | ± 3.5 V @ 50 Ω | ± 3.5 V @ 50 Ω | ± 1 V @ 50 Ω | ± 1 V @ 50 Ω | ± 1 V @ 50 Ω |

Output: 50 Ω, BNC. Offset adjustable by trimpot. Output short-circuit protected. Power supply ± 15 V via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.

Variable Gain 100/200 MHz Voltage Amplifiers ■ Series DHPVA



Model DHPVA-200

| Model | DHPVA-100 | DHPVA-200 |
|-------------------------|-----------------------------|-----------------------------|
| Lower Cut-Off Frequency | DC/10 Hz | DC/10 Hz |
| Upper Cut-Off Frequency | 10/100 MHz, Switchable | 20/200 MHz, Switchable |
| Gain | 10/20/30/40/50/60 dB | 10/20/30/40/50/60 dB |
| Gain Accuracy | ± 0.3 dB | ± 0.3 dB |
| Input Noise Voltage | 2.5 nV/√Hz | 2.5 nV/√Hz |
| Input Voltage Drift | 0.6 μV/°C | 0.6 μV/°C |
| Input/Output | 50 Ω, BNC | 50 Ω, BNC |
| Output Voltage (Power) | ± 1 V (+10 dBm) @ 50 Ω Load | ± 1 V (+10 dBm) @ 50 Ω Load |
| Monitor Output | DC - 100 kHz | DC - 100 kHz |
| Control Interface | TTL/CMOS, Opto-Isolated | TTL/CMOS, Opto-Isolated |

Bandwidth and rise/fall times are independent of gain setting – guaranteed and 100% tested. Offset adjustable by trimpot or external analog control voltage. LED gain setting indication. Output short-circuit protected. Power supply ± 15 V via 3-pin-LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.

- Variable Gain from 10 dB to 60 dB (x 3 to x 1,000)
- Bandwidth DC to 200 MHz for all Gain Settings
- Exceptional Low DC Drift of just 0.6 μV/°C
- Switchable AC/DC Coupling
- Adjustable Offset
- Ideal as Pre-Amplifier for Oscilloscopes, A/D Converters and Transient Recorders

Logarithmic Wideband Voltage Amplifier ■ Model HLVA-100




Model HLVA-100

| Model | HLVA-100 |
|----------------------|---|
| Input Voltage Range | Switchable, ± 20 μV ... ± 200 mV and ± 200 μV ... ± 2 V |
| Dynamic Range | Typ. 60 dB (for Accurate Amplitude Measurement), Max. 80 dB (Signal Detection) |
| Scaling | 12.5 mV/dB, 250 mV/Decade @ 50 Ω Load |
| Linearity | ± 1 dB (for Pulse of Min. 20 ns Pulse Width) |
| Signal Path | DC Coupled, Rectifying |
| Input Noise Voltage | 2 nV/√Hz |
| Input Voltage Drift | 0.6 μV/°C |
| Input/Output | 50 Ω, BNC |
| Rise/Fall Time | 5 ns @ 40 dB Step |
| Output Voltage Range | + 50 ... + 1,075 mV Typ. @ 50 Ω Load (if Output Is Adjusted to 1 V at 100 mV Input) |
| Offset Voltage Range | ± 2.5 mV on Input, ± 500 mV on Output |
| Control Interface | TTL/CMOS, Opto-Isolated |

Offset adjustable by trimpot or external control voltage. Power supply ± 15 V via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.

- Wide Dynamic Range up to 80 dB
- 5 ns Rise/Fall Time @ 40 dB Step
- Switchable Input Range
- DC Coupled, Rectifying Input
- 2 nV/√Hz Input Noise
- Integrated Sample & Hold Baseline Correction
- Typical Applications: LIDAR Systems, Signal Compression, Time-Resolved Pulse and Transient Measurements

High Speed GHz Amplifiers ■ Series HSA-X



Model HSA-X-2-40


- Bandwidth 10 kHz up to 2 GHz
- Fixed Gain of 20 or 40 dB (x 10 or x 100)
- Very Low Input Noise down to 330 pV/√Hz (6.6 pA/√Hz)
- Integrated DC Path for Easy Operation with Fast Photodiodes
- Ultra Compact Package for Use Close to the Signal Source
- Designed for High Speed Photodiodes, APDs, PMTs and as Pre-Amplifier for Fast Oscilloscopes and Transient Recorders

| Model | HSA-X-2-20 | HSA-X-2-40 | HSA-X-1-40 |
|----------------------------|--|------------------------|-------------------------|
| Lower Cut-Off Frequency | 10 kHz | 10 kHz | 10 kHz |
| Upper Cut-Off Frequency | 2 GHz | 2 GHz | 1.1 GHz |
| Rise/Fall Time (10% - 90%) | 180 ps | 180 ps | 320 ps |
| Gain/Transimpedance | 20 dB / 500 V/A | 40 dB / 5,000 V/A | 40 dB / 5,000 V/A |
| Input Noise | 650 pV/√Hz (13 pA/√Hz) | 670 pV/√Hz (13 pA/√Hz) | 330 pV/√Hz (6.6 pA/√Hz) |
| Input VSWR | 1.1 | 1.1 | 1.3 |
| Output Voltage | 1.9 Vp-p (@ 50 Ω Load) | 1.9 Vp-p (@ 50 Ω Load) | 2 Vp-p (@ 50 Ω Load) |
| Input/Output | 50 Ω, SMA | 50 Ω, SMA | 50 Ω, SMA |
| Power Supply | + 15 V, + 85 mA Typ. | + 15 V, + 125 mA Typ. | + 15 V, + 120 mA Typ. |
| Case | 80 x 45 x 25 mm (L x W x H), Weight 100 g (0.23 lbs) | | |

Integrated DC Path. 8-32 and M4 mounting threads. Power supply via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.

High Speed GHz Amplifiers ■ Series HSA-Y

- Bandwidth 10 kHz up to 2 GHz
- Fixed Gain of 20, 40 or 60 dB (x 10, x 100 or x 1,000)
- Very Low Input Noise down to 330 pV/√Hz (6.6 pA/√Hz)
- Two Identical Signal Outputs to Avoid External HF Signal Splitters
- Additional DC Monitor Output for Measuring Low Frequency Background
- Integrated DC Path for Easy Operation with Fast Photodiodes
- Designed for High Speed Photodiodes, APDs, PMTs and as Pre-Amplifier for Fast Oscilloscopes and Transient Recorders




Model HSA-Y-1-60

| Model | HSA-Y-2-20 | HSA-Y-2-40 | HSA-Y-1-40 | HSA-Y-1-60 |
|----------------------------|--|------------------------|-------------------------|-------------------------|
| Lower Cut-Off Frequency | 10 kHz | 10 kHz | 10 kHz | 10 kHz |
| Upper Cut-Off Frequency | 2 GHz | 1.9 GHz | 1.1 GHz | 1.1 GHz |
| Rise/Fall Time (10% - 90%) | 180 ps | 185 ps | 320 ps | 320 ps |
| Gain/Transimpedance | 20 dB / 500 V/A | 40 dB / 5,000 V/A | 40 dB / 5,000 V/A | 60 dB / 50,000 V/A |
| Input Noise | 680 pV/√Hz (14 pA/√Hz) | 650 pV/√Hz (13 pA/√Hz) | 330 pV/√Hz (6.6 pA/√Hz) | 330 pV/√Hz (6.6 pA/√Hz) |
| Input VSWR | 1.15 | 1.2 | 1.45 | 1.4 |
| Input | 50 Ω, SMA | 50 Ω, SMA | 50 Ω, SMA | 50 Ω, SMA |
| Signal Output | Two Identical Signal Outputs, 50 Ω, SMA, Output Voltage 2 Vp-p (@ 50 Ω Load) | | | |
| Monitor Output | Gain 1,000 V/A, Output Voltage: +/- 10 V (@ ≥ 1 MΩ Load), Bandwidth: DC - 100 kHz, BNC | | | |
| Power Supply | ± 15 V, + 190 mA / - 10 mA Typ. | | | |
| Case | 110 x 70 x 25 mm (L x W x H), Weight 180 g (0.41 lbs) | | | |

Integrated DC Path. 8-32 and M4 mounting threads. Power supply via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.

Variable Gain GHz Amplifiers ■ Series DUPVA



Model DUPVA-1-70

- Bandwidth 1 kHz to 1 GHz for all Gain Settings
- Variable Gain from 20 dB to 70 dB (x 10 to x 3,000)
- Very Low Input Noise of 330 pV/√Hz
- Exceptional Gain Flatness of ± 0.15 dB
- Ideal as Broadband Pre-Amplifier for Oscilloscopes and Transient Recorders

| Model | DUPVA-1-60 | DUPVA-1-70 |
|----------------------------|--|----------------------------------|
| Lower Cut-Off Frequency | 1 kHz | 1 kHz |
| Upper Cut-Off Frequency | 1.2 GHz | 1.1 GHz |
| Rise/Fall Time (10% - 90%) | 380 ps | 390 ps |
| Gain | 20/30/40/50/60 dB | 30/40/50/60/70 dB |
| Gain Flatness | ± 0.15 dB | ± 0.15 dB |
| Input Noise | NF 3.0 dB (450 pV/√Hz) | NF 1.9 dB (330 pV/√Hz) |
| Output Power | 10 dBm (-1 dB Compression) | 11 dBm (-1 dB Compression) |
| Output Voltage | 1.7 Vp-p (@ 50 Ω Load) | 1.7 Vp-p (@ 50 Ω Load) |
| Input/Output | 50 Ω, SMA | 50 Ω, SMA |
| Power Supply | ± 15 V, + 350 mA / - 100 mA Typ. | ± 15 V, + 250 mA / - 100 mA Typ. |
| Monitor Output | DC ... 100 kHz, Monitor Output at D-Sub Connector, Gain x1 | |
| Control Interface | 3 Opto-Isolated Digital Inputs, TTL/CMOS Compatible | |
| Case | 165 x 105 x 45 mm (L x W x H), Weight 510 g (1.2 lbs) | |

LED gain setting indication. Output short-circuit protected. Power supply via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.

Miniature Lock-In Amplifier Modules ■ Series LIA-MV-150



Model LIA-MV-150-S

- Working Frequency up to 45 kHz
- Single Ended or True Differential Voltage and Current Input with Sensitivity down to 3 μ V and 3 nA for Full Scale Output
- Phase Shifter 0° to 360°
- Manual and Remote Control Operation
- Ultra Compact Package
- For Laser Stabilization, Spectroscopy, Surface Analysis and Process Control in Scientific and Industrial Applications

| Model | LIA-MV-150-S Standard | LIA-MV-150-D True Differential Input |
|--------------------------|---|--|
| Voltage Input | BNC, Single Ended Instrumentation Amplifier, Noise 12 nV/ \sqrt Hz | LEMO, True Differential Instrumentation Amplifier, Noise 12 nV/ \sqrt Hz |
| Current Input | BNC, Transimpedance Amplifier, Gain 1 kV/A, Noise 13 pA/ \sqrt Hz | LEMO, Transimpedance Amplifier, Gain 1 kV/A, Noise 13 pA/ \sqrt Hz |
| Sensitivity (Full Scale) | Voltage: 3 μ V - 100 mV, Switchable in 1-3-10 Steps / Current: 3 nA - 100 μ A, Switchable in 1-3-10 Steps | |
| Gain Accuracy | \pm 2% for a Sinusoidal Input Signal | |
| Working Frequency | 10 Hz - 45 kHz | |
| Reference Input | BNC, \pm 100 mV to \pm 5 V, Switchable to TTL | |
| Phase | Adjustable 0° - 360° | |
| Demodulator | 55 dB Max. Dynamic Reserve, Square Wave Mixer | |
| Time Constants | 3 ms - 10 s, Switchable in 1-3-10 Steps, Slope Switchable 6 or 12 dB/Octave | |
| Output | BNC, X (In Phase), \pm 10 V Full Scale, Short-Circuit Protected | |
| Control Interface | 16 Opto-Isolated Digital Inputs, TTL/CMOS Compatible, 8 Bit Phase, 4 Bit Time Constant, 4 Bit Sensitivity | |
| Power Supply | \pm 15 V, + 100 mA / - 60 mA Typ. | |
| Case | 170 x 60 x 30 mm (L x W x H), Weight 370 g (0.82 lbs) | |

LED overload and unlocked indication. Power supply via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.

Lock-In Amplifier Modules ■ Series LIA-MV(D)-200

| Model | LIA-MV-200-L Single Phase | LIA-MV-200-H Single Phase | LIA-MVD-200-L Dual Phase | LIA-MVD-200-H Dual Phase |
|--------------------------|--|------------------------------|-----------------------------|-----------------------------|
| Working Frequency | 5 Hz ... 10 kHz | 50 Hz ... 120 kHz | 5 Hz ... 10 kHz | 50 Hz ... 120 kHz |
| Time Constants | 3 ms - 10 s | 300 μ s - 1 s | 3 ms - 10 s | 300 μ s - 1 s |
| Outputs | X = In Phase, Y = Quadrature, R = Magnitude | | | |
| Sensitivity (Full Scale) | Voltage: 3 μ V - 1 V in 1-3-10 Steps / Current: 30 pA - 10 μ A in 1-3-10 Steps | | | |
| Gain Accuracy | \pm 2% for a Sinusoidal Input Signal | | | |
| Voltage Input | BNC, True Differential Instrumentation Amplifier, Noise 12 nV/ \sqrt Hz | | | |
| Current Input | BNC, Transimpedance Amplifier, Gain 100 kV/A, Noise 0.4 pA/ \sqrt Hz | | | |
| Signal Filter | High Pass 0.2 Hz - 1 kHz and Low Pass 100 Hz - 1 MHz, User Selectable | | | |
| Reference Input | BNC, \pm 100 mV to \pm 5 V, Switchable to TTL | | | |
| Phase | Adjustable 0° - 360° | | | |
| Demodulator | 80 dB Max. Dynamic Reserve, Square Wave Mixer | | | |
| Control Interface | 16 Opto-Isolated Digital Inputs, TTL/CMOS Compatible | | | |
| Case | 230 x 105 x 65 mm (L x W x H), Weight 1 kg (2.2 lbs) | | | |

LED overload, unlocked and power indication. Connectors for an optional reference oscillator are inside the module. Power supply \pm 15 V via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.



Model LIA-MVD-200-H

- Single and Dual Phase Models
- Working Frequency 5 Hz up to 120 kHz
- Various Signal Filters and Configuration Options for Versatile Use
- Manual And Remote Control Operation
- Designed as Alternative to Expensive Desktop Lock-Ins for Use in Spectroscopy, Photonics and Laser Measurements

Single-Board Lock-In Amplifiers ■ Series LIA-BV(D)-150



Model LIA-BVD-150-H with Optional Mounting Kit MK-LIA-2

- Low-Cost 19" Plug-In Board Design
- Single and Dual Phase Models
- Working Frequency 5 Hz up to 120 kHz
- Phase Shifter 0° to 360°
- Voltage and Current Input
- Manual and Remote Control Operation
- Ideal for Multi-Channel, OEM and Cost-Sensitive Applications

| Model | LIA-BV-150-L | LIA-BV-150-H | LIA-BVD-150-L | LIA-BVD-150-H |
|---------------------------------|--|-------------------|-----------------|-------------------|
| | Single Phase | Single Phase | Dual Phase | Dual Phase |
| Working Frequency | 5 Hz ... 10 kHz | 50 Hz ... 120 kHz | 5 Hz ... 10 kHz | 50 Hz ... 120 kHz |
| Time Constants | 3 ms – 10 s | 300 μs – 1 s | 3 ms – 10 s | 300 μs – 1 s |
| Outputs | X = In Phase ± 10 V Full Scale, Short-Circuit Protected, Signal Monitor Output X = In Phase, Y = Quadrature, R = Magnitude | | | |
| Sensitivity (Full Scale) | Voltage: 3 μV – 1 V in 1-3-10 Steps / Current: 30 pA – 10 μA in 1-3-10 Steps | | | |
| Gain Accuracy | ± 2% for a Sinusoidal Input Signal | | | |
| Voltage Input | True Differential Instrumentation Amplifier, Noise 12 nV/√Hz | | | |
| Current Input | Transimpedance Amplifier, Gain 100 kV/A, Noise 0.4 pA/√Hz | | | |
| Signal Filter | High Pass 0.2 Hz – 1 kHz and Low Pass 100 Hz – 1 MHz, User Selectable | | | |
| Reference Input | ± 100 mV to ± 5 V, Switchable to TTL | | | |
| Phase | Adjustable 0° - 360° | | | |
| Demodulator | 80 dB Max. Dynamic Reserve, Square Wave Mixer | | | |
| Control Interface | 16 Opto-Isolated Digital Inputs, TTL/CMOS Compatible | | | |
| Dimensions | 160 x 100 x 20 mm (L x W x H), Weight 100 g (0.22 lbs) | | | |

LED overload, unlocked and power indication. Power supply ± 15 V. Connectors for an optional reference oscillator module are on board. Optional mounting kit MK-LIA-2 available. For further information please view the datasheet at www.femto.de.

Multi-Channel Lock-In Amplifier Rack ■ Series SC-LIA-S


| Model | SC-LIA-S |
|--------------------------|---|
| Case Width | 28, 42, 63 or 84 TE |
| Case Height | 3 HE |
| Channel Count | Min. 3, Max. 10 LIA-BV(D)-150 Boards |
| Connections | Signal In-/Outputs: BNC; Control Interface: D-Sub 25 (if Installed) |
| Control Interface | Either Common for All Installed Channels, or Individual for Each Channel, or None |
| Power Supply | Input Voltage: 230 VAC or 110 VAC, Output Voltage: ± 15 VDC, Linear Regulated |
| Case | Depending on the Configuration up to ca. 370 mm x 490 mm x 140 mm (L x W x H), Weight up to ca. 8 kg (17 lbs) |

- 19" Rack System for Series LIA-BV(D)-150 Lock-In Amplifier Boards
- Different Rack Sizes Suitable for 3 up to 10 Channels
- Optional Interface for Remote Control of the Installed LIA Boards
- Integrated Power Supply
- BNC Sockets for Signal In- and Outputs
- Ideal for Applications Requiring Multi-Channel Lock-In Detection



Model SC-LIA-S

Selection Guide ■ Photoreceivers



| Model | Spectral Range | Calibration Wavelength | Bandwidth (-3 dB) | Min. Rise Time (10% - 90%) | Max. Conversion Gain | Min. NEP |
|----------------|-----------------|------------------------|--------------------|----------------------------|--------------------------|------------|
| FWPR-20-SI | 320 ... 1100 nm | - | DC ... 20 Hz | 18 ms | 0.6×10^{12} V/W | 0.7 fW/√Hz |
| FWPR-20-IN | 900 ... 1700 nm | - | DC ... 20 Hz | 18 ms | 9.5×10^{10} V/W | 7.5 fW/√Hz |
| LCA-S-400K-SI | 400 ... 1100 nm | - | DC ... 400 kHz | 1 μs | 6.2×10^6 V/W | 130 fW/√Hz |
| LCA-S-400K-IN | 900 ... 1700 nm | - | DC ... 400 kHz | 1 μs | 9.5×10^6 V/W | 75 fW/√Hz |
| OE-200-SI | 320 ... 1060 nm | 830 nm* | DC ... 500 kHz | 700 ns | 1.0×10^{11} V/W | 10 fW/√Hz |
| OE-200-UV | 190 ... 1000 nm | 830 nm* | DC ... 500 kHz | 700 ns | 1.0×10^{11} V/W | 17 fW/√Hz |
| OE-200-IN1 | 900 ... 1700 nm | 1300 nm* | DC ... 500 kHz | 700 ns | 1.0×10^{11} V/W | 8 fW/√Hz |
| OE-200-IN2 | 900 ... 1700 nm | 1550 nm* | DC ... 500 kHz | 700 ns | 1.0×10^{11} V/W | 7 fW/√Hz |
| HCA-S-200M-SI | 320 ... 1000 nm | - | DC ... 200 MHz | 1.8 ns | 1.1×10^4 V/W | 9.3 pW/√Hz |
| HCA-S-200M-IN | 900 ... 1700 nm | - | DC ... 200 MHz | 1.8 ns | 1.9×10^4 V/W | 5.4 pW/√Hz |
| HCA-S-400M-SI | 320 ... 1000 nm | - | DC ... 400 MHz | 1.0 ns | 2.7×10^3 V/W | 40 pW/√Hz |
| HCA-S-400M-IN | 900 ... 1700 nm | - | DC ... 400 MHz | 1.0 ns | 4.8×10^3 V/W | 24 pW/√Hz |
| HSA-X-S-1G4-SI | 320 ... 1000 nm | - | 10 kHz ... 1.4 GHz | 250 ps | 2.5×10^3 V/W | 26 pW/√Hz |
| HSA-X-S-2G-IN | 900 ... 1700 nm | - | 10 kHz ... 2 GHz | 180 ps | 4.8×10^3 V/W | 14 pW/√Hz |

* Optical calibration available for fiber optic models only.

Femtowatt Photoreceiver ■ Series FWPR-20



Model FWPR-20-IN
Post holder and post not included

- **Ultra Low Noise:**
Min. NEP 0.7 fW/√Hz Allows Direct Detection down to 50 fW
- **Ultra High Gain Amplifier with**
Transimpedance of up to 10¹² V/A
- **For Fluorescence Measurements,**
Spectroscopy, Chromatography,
Electrophoresis and as Replacement
for Photomultiplier Tubes (PMTs) and
Avalanche Photodiodes (APDs)

| Model | FWPR-20-SI | FWPR-20-IN |
|-----------------------|---------------------------------------|---|
| Spectral Range | 320 ... 1100 nm | 900 ... 1700 nm |
| Bandwidth (-3 dB) | DC ... 20 Hz | DC ... 20 Hz |
| Rise Time (10% - 90%) | 18 ms | 18 ms |
| Transimpedance Gain | 1 x 10 ¹² V/A | 1 x 10 ¹¹ V/A |
| Max. Conversion Gain | 0.6 x 10 ¹² V/W (@ 960 nm) | 0.95 x 10 ¹¹ V/W (@ 1550 nm) |
| Min. NEP | 0.7 fW/√Hz (@ 960 nm) | 7.5 fW/√Hz (@ 1550 nm) |
| Saturation Power | 18 pW (@ 960 nm) | 110 pW (@ 1550 nm) |
| Detector | Si, 1.1 x 1.1 mm ² | InGaAs PIN, Ø 0.5 mm |
| Input | Free Space, 25 mm Ø Flange | |
| Output | BNC | |
| Output Voltage Range | ± 10 V @ ≥ 1 MΩ Load | |
| Power Requirements | ± 15 V, ± 15 mA Typ. | |

Threaded M4 and 8-32 holes for mounting on standard posts. 25 mm Ø flange compatible with microbench systems. Offset adjustable by trimpot. Fiber optic input optional. Output short-circuit protected. Power supply via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.

400 kHz Low Noise Photoreceiver ■ Series LCA-S-400K

- **Low Noise:** Min. NEP 75 fW/√Hz Allows Detection down to 1 nW
- **High Gain:** Max. 9.5 x 10⁶ V/W
- **Broad Wavelength Range:** 400 to 1700 nm
- **For Spectroscopy, General Purpose**
Opto-Electronic Measurements and as
Optical Front-End for Oscilloscopes,
A/D Converters and Lock-In Amplifiers

| Model | LCA-S-400K-SI | LCA-S-400K-IN |
|-----------------------|--------------------------------------|---------------------------------------|
| Spectral Range | 400 ... 1100 nm | 900 ... 1700 nm |
| Bandwidth (-3 dB) | DC ... 400 kHz | DC ... 400 kHz |
| Rise Time (10% - 90%) | 1 µs | 1 µs |
| Transimpedance Gain | 1 x 10 ⁷ V/A | 1 x 10 ⁷ V/A |
| Max. Conversion Gain | 6.2 x 10 ⁶ V/W (@ 900 nm) | 9.5 x 10 ⁶ V/W (@ 1550 nm) |
| Min. NEP | 130 fW/√Hz (@ 900 nm) | 75 fW/√Hz (@ 1550 nm) |
| Saturation Power | 1.6 µW (@ 900 nm) | 1 µW (@ 1550 nm) |
| Detector | Si PIN, Ø 2.5 mm | InGaAs PIN, Ø 0.5 mm |
| Input | Free Space, 25 mm Ø Flange | |
| Output | BNC | |
| Output Voltage Range | ± 10 V @ ≥ 1 MΩ Load | |
| Power Requirements | ± 15 V, ± 40 mA Typ. | |

Threaded M4 and 8-32 holes for mounting on standard posts. 25 mm Ø flange compatible with microbench systems. Offset adjustable by trimpot. Fiber optic input and AC-coupling optional. Output short-circuit protected. Power supply via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.



Model LCA-S-400K-SI
Post holder and post not included

500 kHz Variable Gain Photoreceiver ■ Series OE-200



Model OE-200-IN1

- Variable Gain over 9 Decades Allows Detection of fW to mW
- High Speed: Min. Rise Time 700 ns, Max. Bandwidth 500 kHz
- Low Noise: NEP down to 7 fW/√Hz
- Broad Spectral Range: 190 to 1700 nm
- Optical Calibration for Fiber Optic Models
- Applications: All Purpose Lab Photoreceiver, Fiber Alignment Systems, Fast Power Monitoring, Linearity Measurements over 10 Decades, Industrial Control Systems

| Model | OE-200-SI | OE-200-UV | OE-200-IN1 | OE-200-IN2 |
|----------------------------|--|--|--|--|
| Spectral Range | 320 ... 1060 nm | 190 ... 1000 nm | 900 ... 1700 nm | 900 ... 1700 nm |
| Calibration Wavelength | 830 nm | 830 nm | 1300 nm | 1550 nm |
| Max. Bandwidth (-3 dB) | DC ... 500 kHz | DC ... 500 kHz | DC ... 500 kHz | DC ... 500 kHz |
| Min. Rise Time (10% - 90%) | 700 ns | 700 ns | 700 ns | 700 ns |
| Conversion Gain Range | 10 ³ - 10 ¹¹ V/W | 10 ³ - 10 ¹¹ V/W | 10 ³ - 10 ¹¹ V/W | 10 ³ - 10 ¹¹ V/W |
| Min. NEP | 10 fW/√Hz | 17 fW/√Hz | 8 fW/√Hz | 7 fW/√Hz |
| Optical Power Input Range | ca. 100 fW ... 2 mW | ca. 200 fW ... 2 mW | ca. 100 fW ... 2 mW | ca. 100 fW ... 2 mW |
| Detector | Si PIN, ø 1.2 mm | Si PIN, 1.1 x 1.1 mm ² | InGaAs PIN, ø 0.1 mm | InGaAs PIN, ø 0.1 mm |
| Input | Free Space, FC, SMA | Free Space, FC, SMA | Free Space, FC | Free Space, FC |
| Output | BNC | | | |
| Output Voltage Range | ± 10 V @ ≥ 1 MΩ Load | | | |
| Accuracy | ± 1 % Electrical, ± 5 % Electro Optical for Fiber Optic Models, ± 15 % Electro Optical for Free Space Models | | | |
| Lowpass Filter | Switchable to 10 Hz | | | |
| Power Requirements | ± 15 V, + 150 mA / - 100 mA Typ. | | | |
| Control Interface | 5 Opto-Isolated Digital Inputs, TTL/CMOS Compatible, Analog Offset Control Voltage Input | | | |
| Case | 170 x 60 x 45 mm (L x W x H), Weight 320 g (0.74 lbs) | | | |

Offset adjustable by trimpot or external control voltage. LED overload indication. Output short-circuit protected. Power supply via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.

200 MHz High Speed Photoreceiver ■ Series HCA-S-200M

- High Speed: Min. Rise Time 1.8 ns
- Low Noise: Min. NEP 5.4 pW/√Hz Allows Detection down to 1 μW
- High Gain: Max. 1.9 x 10⁴ V/W
- Applications: Spectroscopy, Optical Triggering, Fast Pulse and Transient Measurements

| Model | HCA-S-200M-SI | HCA-S-200M-IN |
|-----------------------|--------------------------------------|---------------------------------------|
| Spectral Range | 320 ... 1000 nm | 900 ... 1700 nm |
| Bandwidth (-3 dB) | DC ... 200 MHz | DC ... 200 MHz |
| Rise Time (10% - 90%) | 1.8 ns | 1.8 ns |
| Transimpedance Gain | 2 x 10 ⁴ V/A | 2 x 10 ⁴ V/A |
| Max. Conversion Gain | 1.1 x 10 ⁴ V/W (@ 800 nm) | 1.9 x 10 ⁴ V/W (@ 1550 nm) |
| Min. NEP | 9.3 pW/√Hz (@ 800 nm) | 5.4 pW/√Hz (@ 1550 nm) |
| Saturation Power | 110 μW (@ 800 nm) | 60 μW (@ 1550 nm) |
| Detector | Si PIN, ø 0.8 mm | InGaAs PIN, ø 0.3 mm |
| Input | Free Space, 25 mm ø Flange | |
| Output | 50 Ω, BNC | |
| Output Voltage Range | ± 1.7 V @ 50 Ω Load | |
| Power Requirements | ± 15 V, ± 60 mA Typ. | |

Threaded M4 and 8-32 holes for mounting on standard posts. 25 mm ø flange compatible with microbench systems. Offset adjustable by trimpot. Fiber optic input optional. Output short-circuit protected. Power supply via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.



Model HCA-S-200M-IN
Post holder and post not included

400 MHz High Speed Photoreceiver ■ Series HCA-S-400M



Model HCA-S-400M-SI
Post holder and post not included

- High Speed: Min. Rise Time 1 ns
- High Gain: Max. 4.8×10^3 V/W
- DC Coupling for Precise Pulse Response
- Applications: Spectroscopy, Optical Triggering, Fast Pulse and Transient Measurements

| Model | HCA-S-400M-SI | HCA-S-400M-IN |
|-----------------------|----------------------------|-----------------------|
| Spectral Range | 320 ... 1000 nm | 900 ... 1700 nm |
| Bandwidth (-3 dB) | DC ... 400 MHz | DC ... 400 MHz |
| Rise Time (10% - 90%) | 1.0 ns | 1.0 ns |
| Transimpedance Gain | 5×10^3 V/A | 5×10^3 V/A |
| Max. Conversion Gain | 2.7×10^3 V/W | 4.8×10^3 V/W |
| Min. NEP | 40 pW/√Hz (@ 800 nm) | 24 pW/√Hz (@ 1550 nm) |
| Saturation Power | 400 μW (@ 800 nm) | 200 μW (@ 1550 nm) |
| Detector | Si PIN, ø 0.8 mm | InGaAs PIN, ø 0.3 mm |
| Input | Free Space, 25 mm ø Flange | |
| Output | 50 Ω, BNC | |
| Output Voltage Range | ± 1.5 V @ 50 Ω Load | |
| Power Requirements | ± 15 V, ± 55 mA Typ. | |

Threaded M4 and 8-32 holes for mounting on standard posts. 25 mm ø flange compatible with microbench systems. Offset adjustable by trimpot. Fiber optic input optional. Output short-circuit protected. Power supply via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.

2 GHz High Speed Photoreceiver ■ Series HSA-X-S

- High Gain: Max. 4.8×10^3 V/W
- High Speed: Min. Rise Time 180 ps
- Low Noise: Min. NEP 14 pW/√Hz Allows Detection down to 10 μW
- Ideal for Fast Optical Measurements

| Model | HSA-X-S-1G4-SI | HSA-X-S-2G-IN |
|-----------------------|----------------------------------|-----------------------------------|
| Spectral Range | 320 ... 1000 nm | 900 ... 1700 nm |
| Bandwidth (-3 dB) | 10 kHz ... 1.4 GHz | 10 kHz ... 2 GHz |
| Rise Time (10% - 90%) | 250 ps | 180 ps |
| Transimpedance Gain | 5×10^3 V/A | 5×10^3 V/A |
| Max. Conversion Gain | 2.5×10^3 V/W (@ 760 nm) | 4.8×10^3 V/W (@ 1550 nm) |
| Min. NEP | 26 pW/√Hz (@ 760 nm) | 14 pW/√Hz (@ 1550 nm) |
| AC Saturation Power | 400 μW (@ 760 nm) | 240 μW (@ 1550 nm) |
| Detector | Si PIN, eff. ø 0.8 mm | InGaAs PIN, eff. ø 0.2 mm |
| Input | Free Space, 25 mm ø Flange | |
| Output | 50 Ω, SMA | |
| Output Voltage Range | 2 Vp-p @ 50 Ω Load | |
| Power Requirements | + 15 V, + 130 mA Typ. | |

Threaded M4 and 8-32 holes for mounting on standard posts. 25 mm ø flange compatible with microbench systems. Model with fiber optic input or DC monitor output optionally available. Output short-circuit protected. Power supply via 3-pin LEMO socket. A mating connector is provided with the device. Optional power supply PS-15 available. For further information please view the datasheet at www.femto.de.



Model HSA-X-S-1G4-SI
Post holder and post not included

Power Supply ■ Series PS-15

- Power Supply Compatible with All FEMTO Amplifiers
- European, US and Asian Version Available
- Short-Circuit Protected
- Linear Regulated and Floating Design for Low Ripple

| Model | PS-15-2-L | PS-15-3-L |
|----------------|--|---------------|
| Input Voltage | 210 - 250 VAC | 100 - 135 VAC |
| Input Plug | Euro Plug, DIN 49464 | UL Plug |
| Output Voltage | ± 15 V, + 400 mA, - 250 mA on LEMO Series 1S, 3-Pin Plug | |
| Ripple | 20 mVp-p Typ. | |



Model PS-15-2-L

USB Control Interface ■ Model LUCI-10



Model LUCI-10


- Digital I/O Interface for USB Remote Control of FEMTO Amplifiers
- Supports Opto-Isolation of Amplifier Signal Path from PC USB Port
- Bus-Powered Operation
- System Driver and Application Software Included
- LabVIEW VI's for Easy Integration in LabVIEW Environment
- Ideal for Remote Control of FEMTO Amplifiers in Automated Systems and Measurement Setups

| Model | LUCI-10 |
|-------------------|--|
| Bus Interface | USB 2.0 Compatible |
| Control Interface | 16 Digital Output Lines, 3 Opto-Isolated Digital Input Lines |
| Connector | USB Typ A, D-Sub, Male (25 Pins) |
| Software | System Driver (.dll) Compatible with C/C++, LabVIEW Library, GUI |

For further information please view the datasheet at www.femto.de.



LabVIEW Library



GUI Frontpanel

In addition to our large variety of standard products we offer professional custom-designed modules and complete solutions for measurement systems. Here are some examples of previous work. Please contact us for details.

■ OEM Solutions with Cost and Size Optimized Design




■ Small-Sized Photodetectors



■ Balanced Photoreceivers



■ Dual Channel Amplifiers



■ Miniature Pick-Up Electronics

This Short Form Catalog has been compiled to provide an overview of the amplifier and photoreceiver products available from FEMTO®.

Detailed datasheets are available for each product listed here. Please visit the datasheet download service at our website www.femto.de.

Specifications are subject to change without notice. Information furnished herein is believed to be accurate and reliable. However, no responsibility is assumed by FEMTO Messtechnik GmbH for its use, nor for any infringement of patents or other rights granted by implication or otherwise under any patent rights of FEMTO Messtechnik GmbH.

Product names mentioned in this product overview may also be trademarks used here for identification purposes only.

© FEMTO Messtechnik GmbH 2010. All rights reserved. FEMTO and the FEMTO logo are trademarks or registered trademarks of FEMTO Messtechnik GmbH in the U.S. and/or other countries.



FEMTO Messtechnik GmbH
Paul-Lincke-Ufer 34 · D-10999 Berlin / Germany
Telephone: +49 (0) 30/446 93 86
Fax: +49 (0) 30/446 93 88
e-mail: info@femto.de · <http://www.femto.de>