

Nan☀Integris

Technical Data Sheet

HiPco[®]

Single-Wall Carbon Nanotubes



光技術をサポートする

株式会社オプトサイエンス


<http://www.optoscience.com>

東京本社 〒160-0014 東京都新宿区内藤町1番地 内藤町ビルディング
TEL:03(3356)1064 FAX:03(3356)3466 E-mail:info@optoscience.com
大阪支店 〒532-0011 大阪市淀川区西中島7-7-2 新大阪ビル西館
TEL:06(6305)2064 FAX:06(6305)1030 E-mail:osk@optoscience.com
名古屋営業所 〒450-0002 名古屋市中村区名駅2-37-21 東海ソフトビル
TEL:052(569)6064 FAX:052(569)8064 E-mail:ngo@optoscience.com


Properties

Individual SWNT Diameter ⁱ	~0.8 – 1.2 nm	
Individual SWNT Length ⁱⁱ	~100 – 1000 nm	
Calculated Molecular Weight ⁱⁱⁱ	~3.4x10 ⁵ – 5.2x10 ⁶ Amu	
Color	Black	
Morphology	Dry powder of nanotubes bundled in ropes	
Maximum Density ^{iv}	1.6 g/cm ³	
Bulk Density ^v	~0.1 g/cm ³	
TGA Residue as Fe ^{vi}		
- Raw	<35 wt%	
- Pure	<15 wt%	
- Super Pure	<5 wt%	
TGA 1 st Derivative Peak Temperature	Raw	~350 – 410°C
	Pure	~470 – 490°C
	Super Pure	~510 – 540°C
TGA Onset Temperature	Raw	~350°C
	Pure	~440°C
	Super Pure	~450°C
Maximum Surface Area ^{vii}	1315 m ² /g	
BET Surface Area	~400 – 1000 m ² /g	
Buckypaper Resistance ^{viii}	~0.2 – 2Ω	
Moisture Content	<5 wt%	


SEM




TEM




Raman



TGA Profile



Particle Size Analysis



- i Diameter distribution measured by Unidym from TEM micrographs. Mean diameter ~1.0 nm.
- ii Measured by Unidym using AFM.
- iii Calculated. Lower limit assumes a SWNT with a diameter of 0.8nm and a length of 100nm. $(0.8\text{nm}/0.245\text{nm})(3.1414)(2 \text{ carbon atoms}) = 20$ carbons around the circumference. For every 0.283nm length there are $4 \times 20 = 80$ carbon atoms. $(100\text{nm}/0.283\text{nm})(80)(12.01) = 339,505$ Amu. Assuming 2 significant digits = 3.4×10^5 . Upper limit assumes a SWNT with a diameter of 1.2nm and a length of 1,000nm. $(1.2\text{nm}/0.245\text{nm})(3.1414)(2 \text{ carbon atoms}) = 31$ carbons around the circumference. For every 0.283nm length there are $4 \times 31 = 124$ carbon atoms. $(1000\text{nm}/0.283\text{nm})(124)(12.01) = 5,262,332$ Amu. Assuming 2 significant digits = 5.2×10^6 .



- iv Calculated assuming single-wall nanotubes of diameter 1.0 nm arranged in crystalline “ropes” or “bundles” (inter-wall spacing 0.3 nm).
- v Value provided is for standard purified SWNTs. Raw and some super pure grades lots will have lower bulk densities. Other product forms may have higher bulk densities.
- vi 800°C in air. The reported figures assume that the residue is present in the product as elemental Fe, and that it is fully converted to Fe_2O_3 during the TGA analysis. Hence, the TGA residual as measured is multiplied by $\text{MW Fe}_2/\text{MW Fe}_2\text{O}_3$ (1/1.43) to express the result as Fe.
- vii Calculated using geometric arguments assuming an isolated tube. SSA for tubes in “ropes” will be less than the stated value. A. Peigney et al., Carbon 39 (2001), 507-514.
- viii In-house, Unidym buckypaper conductivity test (4 point probe).
- ix The following particle size histogram is indicative of the typical tertiary particle size distribution found in bulk powder. It does not represent secondary particles (aggregates of individual tubes also known as ropes or bundles) nor does it represent primary particle sizes (individual carbon nanotubes).